

## **Method Validation Worksheets**

Refer to QA Manual for instructions on using the Method Validation Worksheets. Also, visit the Qualigen, Inc. website at [www.qualigeninc.com](http://www.qualigeninc.com) for “On Q” Training.

**1. Verify Accuracy and Precision – this is performed once per analyzer per assay during installation and training.**

**2. Verify Reportable Ranges (Calibration Verification) – this is performed initially at start up and routinely every 6 months thereafter.**

**3. Identify your Reference Ranges**

**SHBG - Verify Accuracy and Precision.** Run C1 and C2 five times each upon installation.

**C1 (Low Control Range)**

**C2 (High Control Range)**

**1**

**Control C1**

*Write "C1", the control kit lot number, and your initials on the FastPack® IP peel-off label and place it here.*

**Control C2**

*Write "C2", the control kit lot number, and your initials on the FastPack® IP peel-off label and place it here.*

**2**

**Control C1**

*Write "C1", the control kit lot number, and your initials on the FastPack® IP peel-off label and place it here.*

**Control C2**

*Write "C2", the control kit lot number, and your initials on the FastPack® IP peel-off label and place it here.*

**3**

**Control C1**

*Write "C1", the control kit lot number, and your initials on the FastPack® IP peel-off label and place it here.*

**Control C2**

*Write "C2", the control kit lot number, and your initials on the FastPack® IP peel-off label and place it here.*

**4**

**Control C1**

*Write "C1", the control kit lot number, and your initials on the FastPack® IP peel-off label and place it here.*

**Control C2**

*Write "C2", the control kit lot number, and your initials on the FastPack® IP peel-off label and place it here.*

**5**

**Control C1**

*Write "C1", the control kit lot number, and your initials on the FastPack® IP peel-off label and place it here.*

**Control C2**

*Write "C2", the control kit lot number, and your initials on the FastPack® IP peel-off label and place it here.*

## Calculating Accuracy and Precision

From the Control Range Card, enter the Upper and Lower limits as well as the Mean of the control material's acceptable range into the appropriate sections below.

### SHBG - C1

Lower Limit

Mean

Upper Limit

Calculate the Average of all 5 values (1+2+...+5)/5:

AVG = \_\_\_\_\_

Calculate the Standard Deviation of all 5 values:

SD = \_\_\_\_\_

Calculate the Install Precision Specification:

IPS = \_\_\_\_\_

$$IPS = 0.8 \times SD \text{ Spec} \quad \left\{ \begin{array}{l} \text{If AVG} < 20, \text{ SD Spec} = 3 \\ \text{If AVG} \geq 20, \text{ SD Spec} = 0.15 \times \text{AVG} \end{array} \right\}$$

Calculate the Chi-Squared Statistic<sup>1</sup> ( $\chi^2_{\text{stat}}$ ):

$\chi^2_{\text{stat}} =$  \_\_\_\_\_

$$\chi^2_{\text{stat}} = SD^2 \times \left( \frac{2.5}{IPS^2} \right)$$

### SHBG - C2

Lower Limit

Mean

Upper Limit

Calculate the average of all 5 values (1+2+...+5)/5:

AVG = \_\_\_\_\_

Calculate the standard deviation (SD) of all 5 values:

SD = \_\_\_\_\_

Calculate the Install Precision Specification:

IPS = \_\_\_\_\_

$$IPS = 0.8 \times SD \text{ Spec} \quad \left\{ \begin{array}{l} \text{If AVG} < 20, \text{ SD Spec} = 3 \\ \text{If AVG} \geq 20, \text{ SD Spec} = 0.15 \times \text{AVG} \end{array} \right\}$$

Calculate the Chi-Squared Statistic ( $\chi^2_{\text{stat}}$ ):

$\chi^2_{\text{stat}} =$  \_\_\_\_\_

$$\chi^2_{\text{stat}} = SD^2 \times \left( \frac{2.5}{IPS^2} \right)$$

Chi-Squared Statistic ( $\chi^2$ ) statistical analysis is an accepted methodology for precision performance evaluations. Refer to CLSI (Clinical and Laboratory Standards Institute) EP5 approved guideline "Evaluation of Precision Performance of Quantitative Measurement Procedures".

## SHBG Accuracy and Precision

**Analyzer SN:**

Circle Your Response

<b>1</b>	Do all control testing values fall within the acceptable QC range?	Y	N	N/A
<b>2</b>	Is the C1 control $\chi^2_{\text{stat}} < \chi^2_{\text{crit}}?$	Y	N	N/A
<b>3</b>	Is the C2 control $\chi^2_{\text{stat}} < \chi^2_{\text{crit}}?$	Y	N	N/A

The critical value for the Chi-Squared Statistic ( $\chi^2_{\text{crit}}$ ) based on 5 measurements and a 95% confidence level is 7.81.

If you can provide a "Yes" answer in all 3 question categories above, check the box below to accept the manufacturer's claims for accuracy and precision. If the answer in any of the above question categories is No, check the box that you DO NOT accept the manufacturer's claims for accuracy and precision and contact Qualigen System Support.

- Accept the manufacturer's claims for accuracy and precision
- DO NOT accept the manufacturer's claims for accuracy and precision

## SHBG

### Verify Reportable Ranges (Calibration Verification) every 6 months

Verify that the FastPack® IP System is accurate to the limits of the reportable range specified by Qualigen, Inc. by using the FastPack® IP SHBG Method Verification Kit.

<p style="text-align: center;"><b>Low</b></p> <p>Target <input data-bbox="313 577 539 634" type="text"/></p>	<p style="text-align: center;"><i>Write "Low Verifier", the lot number, and your initials on the peel-off FastPack® label and place it here</i></p>
<p style="text-align: center;"><b>Mid</b></p> <p>Range <input data-bbox="302 997 539 1054" type="text"/></p>	<p style="text-align: center;"><i>Write "Mid Verifier", the lot number, and your initials on the peel-off FastPack® label and place it here</i></p>
<p style="text-align: center;"><b>High</b></p> <p>Target <input data-bbox="313 1417 539 1474" type="text"/></p>	<p style="text-align: center;"><i>Write "High Verifier", the lot number, and your initials on the peel-off FastPack® label and place it here</i></p>

The Low Verifier results must be  $< 0.80$  nmol/L and the High Verifier must be  $> 174$  nmol/L to accept the Manufacturer's Reportable Range. If either result is different than the above, record the reportable range based on actual observed values. The mid value must be within the range designated on the range card. If it is not, repeat the test. If any value is out of range after repeating the test, contact Qualigen System Support.

- Accept Manufacturer Reportable Range
- DO NOT Accept. Derived Reportable Range: \_\_\_\_\_

**Identify the Reference Ranges for your Practice**

Each laboratory should determine ranges for their local population. A reference interval study employing serum samples from 237 apparently healthy subjects with no known pre-existing endocrine disorders was performed. The nonparametric 2.5th - 97.5th percentiles (central 95%) were determined for reference partitions as shown below.

Partition	N	Median (nmol/L)	Reference Interval (nmol/L)
Males 18 – 50 years	79	25.6	8.4 – 64.3
Males > 50 years	58	34.1	8.8 – 107.4
Females 12 – 46 years	53	37.8	8.7 – 126.6
Females > 46 years	47	51.8	8.9 – 131.1

Type of Patient	Your Reference (Normal) Ranges
Normal	
Others:	

- OK to begin testing       NOT OK to begin testing

\_\_\_\_\_  
Testing Analyst/Technical Consultant

\_\_\_\_\_  
Date

\_\_\_\_\_  
Laboratory Director

\_\_\_\_\_  
Date