

SHBG

Verify Reportable Ranges (Calibration Verification) every 6 months

Verify that the FastPack® IP System is accurate to the limits of the reportable range specified by Qualigen, Inc. by using the FastPack® IP SHBG Method Verification Kit.

<p style="text-align: center;">Low</p> <p>Target <input data-bbox="310 575 537 632" type="text"/></p>	<p style="text-align: center;"><i>Write "Low Verifier", the lot number, and your initials on the peel-off FastPack® label and place it here</i></p>
<p style="text-align: center;">Mid</p> <p>Range <input data-bbox="298 995 537 1052" type="text"/></p>	<p style="text-align: center;"><i>Write "Mid Verifier", the lot number, and your initials on the peel-off FastPack® label and place it here</i></p>
<p style="text-align: center;">High</p> <p>Target <input data-bbox="310 1415 537 1472" type="text"/></p>	<p style="text-align: center;"><i>Write "High Verifier", the lot number, and your initials on the peel-off FastPack® label and place it here</i></p>

The Low Verifier results must be < 0.80 nmol/L and the High Verifier must be > 174 nmol/L to accept the Manufacturer's Reportable Range. If either result is different than the above, record the reportable range based on actual observed values. The mid value must be within the range designated on the range card. If it is not, repeat the test. If any value is out of range after repeating the test, contact Qualigen System Support.

- Accept Manufacturer Reportable Range
- DO NOT Accept. Derived Reportable Range: _____