## FastPack® Vitamin D Calibrator

For use in conjunction with the FastPack® IP Vitamin D Immunoassay and FastPack® IP System analyzer

**CAUTION: United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.**

**In Canada, use of this product is restricted to laboratories only.**

### INTENDED USE

The FastPack® Vitamin D Calibrator is intended to calibrate the FastPack® IP System when used for the quantitative determination of Vitamin D in human serum and plasma.

### SUMMARY AND PRINCIPLE

Quantitative assay calibration is the process by which a set of samples with known analyte concentrations (i.e. assay calibrators) is tested to measure the response. The mathematical relationship between the measured responses and the known analyte concentrations establishes the calibration curve. The FastPack® IP System analyzer must be calibrated by the user to ensure that it is properly adjusted for the particular lot of FastPacks that are being used. The FastPack® Vitamin D Calibrator is used for this purpose.

### PRODUCT INFORMATION

- The FastPack® Vitamin D Calibrator is included in the FastPack® IP Vitamin D Immunoassay Kit Complete – Cat. No. 25000068.
- Provided ready to use.
- Mix contents by gently inverting before use. Avoid bubble formation.
- Calibrator: 3 mL/vial. Liquid. Contains Tris buffer solution with protein stabilizers. Does not contain any analyte. Calibrator A
- Calibration Card: 1
- Preservative: 0.25% ProClin® 300

### WARNINGS AND PRECAUTIONS

- **For *In Vitro* diagnostic use only.**
- Do not pipette by mouth.
- Do not eat, drink, or smoke in designated work areas.
- Do not mix calibrators from different lots.
- After opening, calibrators are stable for 60 days when stored and handled as directed. Do not use calibrators beyond the expiration date.
- Avoid microbial contamination of reagent when removing aliquots from the bottles.
- Refer to the FastPack® IP System Procedure Manual for calibration procedures.
- Discard unused or expired calibrator material, in stoppered vial, into a Biohazard container.
- The components containing ProClin® 300 are classified per applicable European Economic Community (EEC) Directives as: Irritant (Xi). The following are appropriate Risk (R) and Safety (S) phrases for ProClin® 300:
  - R36/38 Irritant to eyes and skin.
  - R43 May cause sensitization by skin contact.
  - S28 After contact with skin, wash immediately with plenty of soap and water.
  - S37 Wear suitable gloves.

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© 2013 Qualigen, Inc. Tous droits réservés. Qualigen et FastPack sont des marques commerciales ou des marques déposées de Qualigen, Inc.
Qualigen, Inc.