The concentration of testosterone in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the testosterone assay method used. Values obtained with different assay methods should not be used interchangeably.

CAUTION: United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.

INTENDED USE

The FastPack® IP Testo and FastPack® Testo Immunoassay is a paramagnetic particle immunoassay for the in vitro quantitative determination of total testosterone in human serum and plasma. It is intended strictly for in-vitro diagnostic use as an aid in the diagnosis and management of conditions involving excess or deficiency of this androgen. The FastPack® IP Testo and FastPack® Testo Immunoassay is designed for use with the FastPack® IP System and the FastPack® System.

SUMMARY

Testosterone (17 β -hydroxyandrost-4-en-3-one) is a C₁₉ steroid hormone and is the major androgen in males. It is predominantly synthesized in the interstitial Leydig cells of the testis and is regulated by interstitial cell stimulating hormone or luteinizing hormone of the anterior pituitary¹. Testosterone is responsible for the masculinization of the genital tract and the development, and maintenance of male secondary sex characteristics, for example, the accessory sex organs, the prostate, seminal vesicles and growth of facial and pubic hair². Testosterone in females is primarily from three sources: the ovaries, the adrenal gland and the peripheral metabolism of the precursor, androstenedione.

Determination of testosterone levels is an integral part in the assessment of hypogonadal states. Some of the major causes of low testosterone levels in males are hypogonadotropic hypogonadism, testicular failure, hyperprolactinemia, Klinefelter's syndrome, hypopituitarism, estrogen therapy and liver cirrhosis. Elevated levels occur in complete androgen resistance (testicular feminization). In females, some of the major causes of increased testosterone levels are ovarian tumors, polycystic ovaries, adrenal tumors and adrenal hyperplasia.⁷ This assay is not intended for pediatric use.

The FastPack® IP Testo and FastPack® Testo immunoassay is based on a competitive immunoassay principle: testosterone in the sample (control or patient) competes with testosterone covalently coupled to a solid phase paramagnetic particle for an alkaline phosphatase labeled anti-testosterone monoclonal antibody. The amount of labeled-anti-testosterone antibody bound to the magnetic particle is inversely proportional to the concentration of total testosterone in the sample.

TEST PRINCIPAL

The FastPack® IP Testo and FastPack® Testo Immunoassay is a competitive chemiluminescence assay.

-] Primary incubation: Antibody solution (a buffer solution containing a monoclonal C3-testosterone-specific antibody labeled with alkaline phosphatase) [100 μ L] reacts with testosterone from the patients sample, control, or calibrator [100 μ L].
-] Secondary incubation: The reaction mixture is added to paramagnetic particle with covalently coupled testosterone. During this incubation, the testosterone-coated beads compete with the sample testosterone.
-] Removal of unbound materials: The paramagnetic particles are washed with wash buffer [0.2 mL/wash] to remove unbound materials.
-] Substrate addition and detection: Chemiluminogenic substrate [140 μ L] is added to the solid-phase bound complex and results in "glow" chemiluminescence, which is measured using the FastPack® IP System and FastPack® System analyzer.
-] The amount of bound labeled-antibody is inversely proportional to the concentration of testosterone in the sample.

REAGENTS - Content and Concentration

Each FastPack® IP carton contains:

-] 30 FastPack® IP

Each FastPack® carton contains:

-] 50 FastPacks

Each FastPack® Contains:

-] Paramagnetic Particles, 150 µL
Paramagnetic particles with covalently coupled testosterone in buffer containing 0.1% sodium azide as a preservative.
-] Testosterone Antibody Solution, 100 µL
Antibody solution containing mouse monoclonal antibody labeled with alkaline phosphatase in a protein matrix containing 0.03% Proclin® 150 as a preservative.
-] Wash Buffer, 2.0 mL
Tris buffer containing surfactants.
-] Substrate, 140 µL
ImmuGlow™: Indoxyl-3-phosphate and lucigenin in buffer containing preservatives.

Materials required but not provided

-] FastPack® IP System or FastPack® System
-] FastPack® Testo Calibrator Kit - Cat. No. 25000014
-] FastPack® Control Kit – Cat. No. 25000056(US) or Cat. No. 25000065(International)

REAGENT WARNING AND PRECAUTIONS

-] For *in vitro* diagnostic use only.
-] Do not pipette by mouth.
-] Do not eat, drink or smoke in designated work areas.
-] Wash hands thoroughly after handling specimen.
-] HAMA Interference: some individuals have antibodies to mouse protein (HAMA), which can cause interference in immunoassays that employ antibodies derived from mice⁵.
-] FastPack® IP reagents and FastPack® reagents are stable until the expiration date on the label when stored and handled as directed. Do not use FastPack® IP reagents and FastPack® reagents beyond expiration date.
-] Discard used FastPacks into a biohazard container.
-] The components containing sodium azide are classified per applicable European Economic Community (EEC) Directives as: Harmful (Xn). The following are appropriate Risk (R) and Safety (S) phrases:
 - R22 Harmful if swallowed.
 - R32 Contact with acids liberates very toxic gas.
 - S2 Keep out of reach of children.
 - S13 Keep away from food, drink and animal feeding stuffs.
 - S36 Wear suitable protective clothing.
 - S46 If swallowed, seek medical advice immediately and show this container or label.

REAGENT STORAGE INSTRUCTIONS

Store at 2 - 8 °C. Protect from light.

SPECIMEN COLLECTION/PREPARATION

1. Serum or lithium heparin or K2 EDTA plasma samples can be used for the FastPack® IP Testo and FastPack® Testo Immunoassay.
2. The National Committee for Clinical Laboratory Standards (NCCLS) provides recommendations for handling, processing, and storing blood.^{3,4}
3. Collect all blood samples observing routine precautions for venipuncture.
4. It is not required that patients fast prior to blood collection.
5. For serum samples:
 -] Collect samples in serum blood collection tubes.
 -] Ensure that complete clot formation has occurred before centrifugation. This takes approximately 30 minutes. Some samples may exhibit increased clotting time, especially those from patients receiving anticoagulant or thrombolytic therapy.
 -] Serum should be centrifuged and separated from the clot within 3 hours from time of collection.
 -] Remove serum from the cells prior to storage at 2-8 °C.
 -] If not tested within 24 hours, the sample should be frozen at –20 °C or colder.
6. For plasma samples:
 -] Collect samples in Lithium Heparin or K2 EDTA tubes.
 -] Plasma should be centrifuged and separated within 3 hours from time of collection.
 -] Remove plasma from the cells prior to storage at 2-8 °C.
 -] If not tested within 24 hours, the sample should be frozen at –20 °C or colder.

7. Do not freeze samples (-20 °C) for more than two months.
8. Frozen samples should be thawed completely and mixed by gentle inversion prior to use.
9. Samples should be free of Fibrin, Red Blood Cells, or other particulate material for optimal results.
10. Samples showing turbidity and/or particulate matter should be centrifuged prior to use.
11. Ensure the samples are free of bubbles.
12. Human samples should be handled in accordance with the OSHA standard on Bloodborne Pathogens.¹⁹

ASSAY PROCEDURE

See the FastPack® IP System Procedure Manual for detailed instructions for running the FastPack® IP System and FastPack® System assays.

INSTRUMENTATION

FastPack® IP System or FastPack® System

DETAILS OF CALIBRATION

During the FastPack® IP and FastPack® production process, Qualigen generates a master standard curve and places this information in the barcode of each FastPack® IP and FastPack® label, where it can be read by the FastPack® IP System or FastPack® System analyzer during the testing sequence. The analyzer must be calibrated by the user to ensure that it is properly adjusted for the particular lot of FastPacks that is being used. Separate calibrations must be run for each type of test, i.e. Free PSA, Total PSA or Testosterone. The frequency of calibration varies for each test type. For the FastPack® IP Testo and FastPack® Testo Immunoassay, the FastPack® IP System and FastPack® System analyzer must be calibrated once every 14 days or whenever a new lot of Testo FastPacks are to be used.

Whenever the user performs an initial calibration for a particular lot of FastPacks or uses a new lot of calibrator, 2 FastPacks must be run for calibration (duplicates). Whenever recalibration is performed with the same lot of FastPacks and calibrator, 2 FastPacks must be run for calibration. See the FastPack® IP System Procedure Manual for "Running a Calibration".

Use FastPack® Testo Calibrator Kit – Cat. No. 25000014

RESULTS

The FastPack® Analyzer uses the information from the barcode to construct a lookup table of x,y values that represent the standard curve and estimates the concentration of unknown samples by linear interpolation.

QUALITY CONTROL

Quality control materials simulate real specimens and are essential for monitoring the system performance of assays. Good Laboratory Practices (GLP) include the use of control specimens to ensure that all reagents and protocols are performing properly. Commercial control products (such as Microgenics Liquimmune® Liquid Assayed Immunoassay Control, level 1 and 2) are available which contain testosterone. See the FastPack® IP System Procedure Manual for "Control Testing".

Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.

Controls available: FastPack® Control Kit – Cat. No. 25000056(US) or Cat. No. 25000065(International)

LIMITATION OF PROCEDURE

-] Samples can be measured within the reportable range of the analytical sensitivity and the upper end of the calibration range, 1600 ng/dL.
-] Samples >1600 ng/dL should be run using another method. Dilution of out of range results is not recommended.
-] Samples from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such samples may show either falsely elevated or depressed values when tested with assay kits employing mouse monoclonal antibodies.
-] Heterophilic antibodies in a sample have the potential to cause interference in immunoassay systems. Infrequently, testosterone levels may appear elevated due to heterophilic antibodies present in the patient's serum or to nonspecific protein binding. If the testosterone level is inconsistent with clinical evidence, additional testosterone testing is suggested to confirm the result.
-] For diagnostic purposes, the Fastpack® IP Testo and FastPack® Testo Immunoassay should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.
-] The FastPack® IP Testo and Fastpack® Testo Immunoassay is not for pediatric use.

EXPECTED RANGE

Serum samples were obtained from 125 random males and 173 random females ranging in age from 18 to 94 years. Samples were obtained from normal, healthy blood donors without any clinically abnormal indications. Testosterone levels were determined using the FastPack® IP Testo and Fastpack® Testo Immunoassay in conjunction with the FastPack® IP System analyzer in order to establish the testosterone concentration in the normal population. Results are in ng/dL.

Test Subjects	N	PERCENTILES				
		2.5th	5th	Median	95th	97.5th
Males 20-49 years	98	214.1	230.6	550.5	939.1	1009.4
Males > 50 years	27	216.2	226.8	439.0	916.0	1004.1
Females 20-49 years	157	<23.0	<23.0	32.0	277.1	318.2
Females > 50 years	16	<23.0	<23.0	29.5	142.8	157.8

Note: The expected range reflects the donor population of this study group. Each laboratory should determine their own reference range appropriate for their population.

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The reproducibility of the testosterone assay was measured by 2 patient serum samples in duplicate (n=60, for each sample) over ten days using three analyzers and two lots of reagents. Precision is expressed either as a Standard Deviation (low level), or Coefficient of Variation (high sample). At the lower level the testosterone value is low thus the CV has little meaning.

Sample	Average (ng/dL)	Between Run	Between Analyzer	Between Reagent Lot	Total
Low Sample	14.0	SD = 6.0	SD = 3.0	SD = 1.0	SD = 7.0
High Sample	712.0	%CV = 10.3	%CV = 3.8	%CV = 5.0	%CV = 10.9

Also, a FastPack® IP Testo and FastPack® Testo immunoassay precision study was conducted per the CLSI EP5-A standard. Three levels of each of the sample types (Lithium Heparin and K2 EDTA plasma) were tested along with 1 level of control. One reagent lot was used for testing. Each test sample was tested 4 reps per day for 20 days.

FastPack® IP Testo and FastPack® Testo EP5-A precision results in ng/dL.

Sample	Average	Within-Day	Total
		%CV	%CV
Control	550.4	7.2	7.3
K2 EDTA Low Pool	72.6	8.6	12.4
K2 EDTA Mid Pool	504.1	7.0	8.4
k2 EDTA High Pool	1092.5	4.4	4.7
Lithium Heparin Low Pool	83.9	8.0	10.5
Lithium Heparin Mid Pool	504.3	7.7	8.8
Lithium Heparin High Pool	1071.0	3.6	3.7

MEASURING RANGE

Spike Recovery:

A patient sample was spiked with testosterone from a human based stock material. The stock concentration was 5400 ng/dL. All samples recover within the $\pm 20\%$ acceptance criteria.

Sample Number	Added Testosterone (ng/dL)	Expected Value Testosterone (ng/dL)	Observed Value Testosterone (ng/dL)	Recovery (%)
1	0	398.0	398.0	100
2	100.0	498.0	533.0	107
3	251.0	649.0	714.0	110
4	501.0	899.0	927.0	103
5	752.0	1150.0	1090.0	94.8
6	1002.0	1400.0	1360.0	97.1

The measuring range was evaluated following the CLSI EP6-A guideline. For Testosterone in serum, Lithium Heparin plasma, and K2 EDTA plasma as tested by the FastPack Testosterone assay, the measuring range is from 23.0 ng/dL (Limit of Quantitation) to 1600 ng/dL. A recovery study was performed to confirm the measuring range. Study results are summarized in the tables below. Diluting samples above 1,600 ng/dL is not recommended.

K2 EDTA

<u>Dilution</u>	<u>Recovery (ng/dL)</u>	<u>Target (ng/dL)</u>	<u>%Bias</u>
1.00	1639.4	1639.4	0.0
0.75	1344.1	1232.8	9.0
0.50	945.7	826.2	14.5
0.25	468.8	419.5	11.7

Lithium Heparin

<u>Dilution</u>	<u>Recovery (ng/dL)</u>	<u>Target (ng/dL)</u>	<u>%Bias</u>
1.00	1686.6	1686.6	0.0
0.75	1342.5	1267.6	5.9
0.50	990.6	848.7	16.7
0.25	517.4	429.7	20.4

Serum

<u>Dilution</u>	<u>Recovery (ng/dL)</u>	<u>Target (ng/dL)</u>	<u>%Bias</u>
1.00	1610.9	1610.9	0.0
0.75	1276.9	1208.2	5.7
0.50	917.9	805.5	13.9
0.25	455.5	402.8	13.1

Method Comparison

Clinical samples were used to compare the serum values obtained using the FastPack® IP Testo and FastPack® Testo method and the serum value obtained using the DPC Coat-a-Count RIA method. The values were evaluated for agreement using Deming regression analysis. Spearman rank correlation was used to test association between the values.

n	Range of Observation (ng/dL)	Intercept (ng/dL)	Slope	R
135	24.0 – 1587.2	6.5	1.054	0.914

MATRIX COMPARISONS: SERUM VS. PLASMA

Lithium Heparin Plasma vs. Serum

Blood collections were obtained from 60 healthy male and female volunteers between the ages of 21 and 56 (samples from 11 of these volunteers were spiked with varying amounts of testosterone to fill out the assay range), and these specimens were processed to heparin plasma and serum samples in parallel.

n	Range of Observation (ng/dL)	Intercept (ng/dL)	Slope	R ²	Sy x
60	24.0 – 1440.0 (Heparin Plasma)				
60	28.0 – 1420.0 (Serum)	-8.0	0.992	0.987	43.9

K2 EDTA Plasma vs. Serum

Blood collections were obtained from 67 healthy male and female volunteers between the ages of 18 and 59 (samples from 11 of these volunteers were spiked with varying amounts of testosterone to fill out the assay range), and these specimens were processed to K2 EDTA plasma and serum samples in parallel.

n	Range of Observation (ng/dL)	Intercept (ng/dL)	Slope	R ²	Sy/x
67	27.5 – 1425.0 (K2 EDTA Plasma)				
67	26.0 – 1420.0 (Serum)	-10.9	0.998	0.996	26.4

INTERFERING SUBSTANCES

Two concentrations of bilirubin, hemoglobin and triglycerides (using Intralipid®) were added to serum, Lithium Heparin plasma, and K2 EDTA plasma samples containing concentrations of <100 ng/dL, approximately 500 ng/dL and approximately 1,000 ng/dL testosterone. The value obtained for the sample with each interfering substance was compared to the value obtained without the interfering substance. Interference was defined as ± 15% recovery for the lowest concentration of interferent and testosterone concentration across all three sample matrices. These compounds did not show interference at the levels indicated below. Triglyceride demonstrated interference for all concentrations of interferent tested at <100 ng/dL of testosterone. Based on these studies, no visibly hemolyzed or cloudy (lipemic) samples should be used with this assay. Lipemia may be removed by ultracentrifugation, if available, and analyzed using the supernatant.

Test Compound	Test Concentration
Bilirubin	10 mg/dL
Hemoglobin	250 mg/dL
Triglycerides	See comment above

Interfering substances were added to control samples containing known amounts of Testo. The value obtained for the control samples with each interfering substance was compared to the value obtained for the control samples without the interfering substance. These compounds did not show interference at the levels indicated.

Test Compound	Test Concentration
d-Biotin	1000 ng/mL

Limit of Blank (LOB), Limit of Detection (LOD), and Limit of Quantitation (LOQ)

The limit of blank (LOB, the highest measurement likely to be observed for a blank sample), limit of detection (LOD, the lowest amount of analyte in a sample that can be distinguished (with 95% confidence) from the blank sample, and limit of quantitation (LOQ, the lowest amount of analyte in a sample that can be reliably detected and at which the total error meets the pre-specified requirement for accuracy) were determined for serum, Lithium Heparin plasma, and K2 EDTA plasma according to CLSI EP17-A. In this study, the Limits of Blank is 3.8 ng/dL, the LOD is 14.5 ng/dL Testosterone, and the LOQ is 23.0 ng/dL Testosterone.

Analytical Specificity

For the monoclonal antibody used, the following cross-reactivity contained in the table below were found. What is tabulated are the maximum % cross-reactivity observed in the study. Twenty three compounds were tested. Serum, Lithium Heparin plasma, and K2 EDTA plasma samples containing concentrations of < 100 ng/dL, approximately 500 ng/dL and approximately 1,000 ng/dL testosterone were spiked with two levels of the cross reactants. % cross-reactivity is defined as 100 x (Testo Recovery – Baseline Testo)/Cross Reactant concentration. Androstenediol and 19-Norethisterone Acetate generated the highest cross-reactivity in this study though although the magnitude was limited to <4%.

Cross-Reactant	Levels Tested (ng/dL)	% Cross-reactivity
5- α -DHT	2500, 5000	2.75
Norethindrone	5000, 10000	1.72
Androstenediol	5000, 10000	3.26
19-Norethisterone Acetate	5000, 10000	3.26
Oxymetholone	5000, 10000	2.37
11-Deoxycortisol	50000, 100000	0.08
19-Nortestosterone	50000, 100000	1.82
Corticosterone	250000, 500000	0.02
Estriol	5000, 10000	0.52
Danazol	50000, 100000	0.02
Estradiol	50000, 100000	0.27
Androstenedione	50000, 100000	1.42
Progesterone	50000, 100000	0.09
Cortisone	50000, 100000	0.24
17-alpha-Hydroxyprogesterone	5000, 10000	1.72
Estrone	25000, 50000	0.22
DHEA Sulfate	50000, 100000	0.08
Ethinylestradiol	5000, 10000	1.56
Methyltestosterone	5000, 10000	1.08
Norgestrel	5000, 10000	1.71
Androsterone	10000000, 50000000	0.02
Cortisol	400000, 800000	0.00
DHEA	500000, 1000000	0.02

REFERENCES

- ¹ Wilson, JD, George, FW, Griffin, JE: The hormonal control of sexual development. *Science*, 211: 1278 – 1284, 1981.
- ² Imperato-Mcginley, J, Guerrero, T, Peterson, RE: Steroid 5 a-reductase deficiency in man: An inherited form of male pseudohermaphroditism. *Science*, 186:1213 – 1215, 1974.
- ³ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition: H3-A5: 23(32) 2003, National Committee for Clinical Laboratory Standards (NCCLS)
- ⁴ Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- ⁵ Schroff, RJ, Foon, KA, et.al.: Human anti-mouse immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res*, 45:879 – 885, 1985
- ⁶ US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. *Federal Register* 1991; 56(235): 64175-82.
- ⁷ Gronowski AM, Landau-Levine M. Reproductive Endocrine Function in Tietz Textbook of Clinical Chemistry. Ed. Edward R. Ashwood, and Carl A. Burtis. 3rd ed. N.p.:W. B. Saunders Company, 1999. 1603-1636.

Qualigen, Inc.
Carlsbad, CA 92011 USA
Technical Support
(760) 918-9165
(877) 709-2169

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© 2000 Qualigen, Inc. All rights reserved. Qualigen and FastPack are trademarks or registered trademarks of Qualigen, Inc. All other trademarks are the property of their respective owners.

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REFERENZEN

- ¹ Wilson, JD, George, FW, Griffin, JE: The hormonal control of sexual development. Science, 211: 1278 – 1284, 1981.
- ² Imperato-Mcginley, J, Guerrero, T, Peterson, RE: Steroid 5 a-reductase deficiency in man: An inherited form of male pseudohermaphroditism. Science, 186:1213 – 1215, 1974.
- ³ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition: H3-A5: 23(32) 2003, National Committee for Clinical Laboratory Standards (NCCLS)
- ⁴ Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- ⁵ Schroff, RJ, Foon, KA, et.al.: Human anti-mouse immunoglobulin responses in patients receiving monoclonal antibody therapy. Cancer Res, 45:879 – 885, 1985
- ⁶ US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. Federal Register 1991; 56(235): 64175-82.
- ⁷ Gronowski AM, Landau-Levine M. Reproductive Endocrine Function in Tietz Textbook of Clinical Chemistry. Ed. Edward R. Ashwood, and Carl A. Burtis. 3rd ed. N.p.:W. B. Saunders Company, 1999. 1603-1636.

Qualigen, Inc.
Carlsbad, CA 92011
USA
Technischer Support
(760) 918-9165
(877) 709-2169

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