For the Quantitative Measurement of Alpha Glutathione S-Transferase ( $\alpha$ GST) in Human serum and plasma

### INTENDED USE

The FastPack® IP  $\alpha$ GST Immunoassay is a paramagnetic particle immunoassay for the *in-vitro* quantitative determination of  $\alpha$ GST in human serum and plasma. Measurement of  $\alpha$ GST can be used as an aid in the assessment of liver damage. The FastPack®  $\alpha$ GST IP Immunoassay is designed for use with the FastPack® IP System.

### SUMMARY

In the liver, alpha glutathione S-transferase is located in the hepatocytes whereas pi GST ( $\pi$ GST) is confined to the intrahepatic bile duct cells. This heterogenous GST subclass distribution suggests that the isoenzymes have unique *in vivo* functions in different hepatic regions and that the detection of GST subclass levels in biological fluids would be of significant use in monitoring the integrity of specific hepatic regions.

Currently, liver injury is studied by the measurement of liver enzymes such as alanine aminotransferase (ALT) and aspartate aminotransferase (AST). A disadvantage of these markers is that they are not distributed uniformly throughout the liver, the periportal concentration being greater than the centrilobular. In contrast,  $\alpha$ GST has been found to be equally distributed in both the centrilobular and periportal regions. Since the centrilobular hepatocytes are very susceptible to damage in a variety of clinical conditions including Allograft Rejection, Viral Hepatitis, Chronic Active Hepatitis and Hepatotoxicity,  $\alpha$ GST would be a more sensitive indicator of hepatic status in these and other, clinical situations. Therefore,  $\alpha$ GST quantitation can be used to study the hepatocellular status of individuals at risk for hepatic damage.<sup>1-7</sup>

### TEST PRINCIPLE

The FastPack®  $\alpha$ GST IP Immunoassay is a chemiluminescence assay based on the “sandwich” principle.

- Primary incubation: Antibody solution (mixture of a biotinylated monoclonal  $\alpha$ GST-specific antibody and a monoclonal  $\alpha$ GST antibody labeled with alkaline phosphatase) [100  $\mu$ L] reacts with  $\alpha$ GST from the patient's sample, control or calibrator [100  $\mu$ L].
- Secondary incubation: Streptavidin-coated paramagnetic particles [150  $\mu$ L] are combined with the reaction mixture. During this incubation, the sandwich complex is bound to the solid-phase via the interaction of biotin and streptavidin.
- Removal of unbound materials: The paramagnetic particles are repeatedly washed with wash buffer [0.2 mL/wash] to remove unbound materials.
- Substrate addition and detection: Chemiluminogenic substrate [140  $\mu$ L] is added to the solid-phase bound complex and results in “glow” chemiluminescence, which is measured using the FastPack® IP System.
- The amount of bound labeled-antibody is directly proportional to the concentration of  $\alpha$ GST in the sample.

### REAGENTS – Content and Concentration

Each FastPack® IP carton contains:

- 30 FastPack® IPs

Each FastPack® IP Contains:

- Paramagnetic Particles, 150  $\mu$ L  
Streptavidin-coated paramagnetic particles in buffer containing 0.1% ProClin® 300 as a preservative.
- $\alpha$ GST Antibody Solution, 100  $\mu$ L  
Antibody solution containing a mixture of a biotinylated mouse monoclonal anti- $\alpha$ GST antibody and a second mouse monoclonal anti- $\alpha$ GST antibody labeled with alkaline phosphatase in a protein matrix containing ProClin® 150 as preservative.
- Wash Buffer, 2.0 mL  
Tris buffer containing surfactants.
- Substrate, 140  $\mu$ L  
ImmuGlow™ Plus: Indoxyl-3-phosphate and lucigenin in buffer containing preservatives.

Materials required but not provided

- FastPack® IP System
- FastPack®  $\alpha$ GST Calibrator Kit – Cat. No. 25000048
- FastPack®  $\alpha$ GST Control Kit – Cat. No. 25000049
- FastPack®  $\alpha$ GST Sample Diluent Kit – Cat. No. 25000050

## WARNINGS AND PRECAUTIONS

- For *in-vitro* diagnostic use only.
- Do not pipette by mouth.
- Do not eat, drink or smoke in designated work areas.
- Wash hands thoroughly after handling specimen.
- HAMA Interference: some individuals have antibodies to mouse protein (HAMA), which can cause interference in immunoassays that employ antibodies derived from mice<sup>9</sup>.
- FastPack<sup>®</sup> IP reagents are stable until the expiration date on the label when stored and handled as directed. Do not use FastPack<sup>®</sup> IP reagents beyond the expiration date.
- Discard used FastPacks into a Biohazard container.
- The components containing ProClin<sup>®</sup> are classified per applicable European Economic Community (EEC) Directives as: Irritant (Xi). The following are appropriate Risk (R) and Safety (S) phrases for ProClin<sup>®</sup>:

R36/38 Irritating to eyes and skin  
R43 May cause sensitization by skin contact  
S24/25 Avoid contact with skin and eyes  
S36/37 Wear suitable protective clothing and gloves  
S60 This material and/or its container must be disposed of as hazardous waste

## STORAGE INSTRUCTIONS

Store at 2 – 8 °C.

## SPECIMEN COLLECTION/PREPARATION

1. Serum, EDTA or lithium-heparin plasma samples can be used for the FastPack<sup>®</sup>  $\alpha$ GST IP Immunoassay.
2. The Clinical and Laboratory Standards Institute provides the following recommendations for handling, processing and storing blood.<sup>9,10</sup>
  - A. Collect all blood samples observing routine precautions for venipuncture.
  - B. For plasma samples:
    - Collect samples in an EDTA (lavender top) or heparinized (green top) tube.
    - Mix the tube immediately after collection by gently inverting it several times.
    - Plasma should be separated from the cells by centrifugation within 3 hours from time of collection and stored at 2 – 8 °C. Transfer the plasma from the original tube for storage.
    - If not tested within 24 hours, the sample should be frozen at –20 °C or colder.
  - C. Samples should be free of red blood cells, or other particulate material for optimal results.
  - D. Samples showing turbidity and particulate matter should be centrifuged prior to use.
  - E. Ensure the samples are free of bubbles.

## ASSAY PROCEDURE

See the FastPack<sup>®</sup> IP System Procedure Manual for detailed instructions for running the FastPack<sup>®</sup> IP assays.

## INSTRUMENTATION

FastPack<sup>®</sup> IP System

## DETAILS OF CALIBRATION

During the FastPack<sup>®</sup> IP production process, Qualigen generates a master standard curve and places this information in the barcode of each FastPack<sup>®</sup> IP label, where it can be read by the FastPack<sup>®</sup> IP System analyzer during the testing sequence. The FastPack<sup>®</sup> IP System analyzer must be calibrated by the end user so that it is properly adjusted for the particular lot of FastPacks that are being used. Separate calibrations must be run for each type of test, i.e. Total PSA, Testosterone or  $\alpha$ GST. The frequency of calibration varies for each test type. For the FastPack<sup>®</sup> IP  $\alpha$ GST Immunoassay, the FastPack<sup>®</sup> IP System analyzer must be calibrated once every 30 days or whenever a new lot of  $\alpha$ GST FastPacks are to be used.

Whenever the user performs a calibration for a particular lot of FastPacks or uses a new lot of calibrator, 2 FastPacks must be run for calibration (duplicates). When the calibration expires (30 days after initial calibration) 2 FastPacks must be run for calibration. See FastPack<sup>®</sup> IP System Procedure Manual for "Running a Calibration".

Use FastPack<sup>®</sup>  $\alpha$ GST Calibrator Kit – Cat. No. 25000048

## RESULTS

The FastPack<sup>®</sup> IP System analyzer uses the information from the barcode to construct a lookup table of x,y values that represent the standard curve and estimates the concentration of unknown samples by linear interpolation.

## QUALITY CONTROL

Quality control materials are essential for monitoring the system performance of assays. Good Laboratory Practices (GLP) includes the use of control specimens to ensure that all reagents and protocols are performing properly. See FastPack<sup>®</sup> IP System Procedure Manual for "Control Testing". At least two levels of quality control materials should be used. Users should follow the appropriate national, state and local guidelines concerning the running of external quality controls.

Controls available: FastPack<sup>®</sup>  $\alpha$ GST Control Kit – Cat. No. 25000049


## LIMITATION OF PROCEDURE

- Plasma samples to be collected using EDTA or lithium-heparin as the anticoagulant.
- Specimens can be measured within the reportable range of the limit of quantitation (0.5 ng/mL) and the upper end of the calibration range of 200 ng/mL.
- Samples >200 ng/mL should be reported as such or re-run after dilution. Samples may be diluted using the  $\alpha$ GST Sample Diluent until values within range of the assay are obtained.
- The FastPack<sup>®</sup> IP  $\alpha$ GST Immunoassay does not show a high-dose hook effect up to 5,000 ng/mL.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits employing mouse monoclonal antibodies.
- Heterophilic antibodies in a sample have the potential to cause interference in immunoassay systems. Infrequently,  $\alpha$ GST levels may appear depressed due to heterophilic antibodies present in the patient's sample or to nonspecific protein binding. If the  $\alpha$ GST level is inconsistent with clinical evidence, additional  $\alpha$ GST testing is suggested to confirm the result.
- For diagnostic purposes, the FastPack<sup>®</sup> IP  $\alpha$ GST Immunoassay should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.
- Performance of this assay has not been established with neonatal specimens.

## EXPECTED VALUES

Samples were obtained from 167 individuals. Samples were obtained from apparently healthy blood donors without any clinically abnormal indications.  $\alpha$ GST levels were determined using the FastPack<sup>®</sup> IP  $\alpha$ GST Immunoassay in conjunction with the FastPack<sup>®</sup> IP System in order to establish the  $\alpha$ GST concentration in the normal population. The reference interval (5<sup>th</sup> to 95<sup>th</sup> percentiles) for the FastPack<sup>®</sup> IP  $\alpha$ GST Immunoassay is 0.24 – 11.4 ng/mL. The reference intervals reflect the donor population of this study group. Each laboratory should determine their own reference range appropriate for their population. The distribution of the values is presented in the histogram below.

Sample Type	Reference Group	n	Median Age (years)	Age Range	Median Concentration (ng/mL)	95% Reference Interval (ng/mL)
Serum	Males and Females	167	35	18 - 55	1.4	0.24-11.4

## SPECIFIC PERFORMANCE CHARACTERISTICS

### Precision

Three control serum pools spiked with varying amounts of  $\alpha$ GST were assayed using three lots of FastPack<sup>®</sup> IP  $\alpha$ GST reagents and four FastPack<sup>®</sup> IP instruments. Each control was tested in duplicate determinations per day over a period of 10 days for a total of 80 replicate determinations per control. The data was analyzed using the analysis of variance (ANOVA) technique and within-assay and total imprecision were determined. Mean, standard deviation (SD), and percent coefficient of variation (%CV) were calculated for each control solution.

#### Within-Assay Imprecision

	Control 1	Control 2	Control 3
Mean ng/mL $\alpha$ GST	7.3	37.6	180.5
SD	0.34	1.5	10.4
%CV	4.7	4.0	5.8

#### Total Imprecision

	Control 1	Control 2	Control 3
Mean mIU/mL $\alpha$ GST	7.6	37.6	180.5
SD	0.80	3.0	25.4
%CV	11.1	7.9	14.1

### Range of Linearity.

For  $\alpha$ GST tested by FastPack<sup>®</sup> IP  $\alpha$ GST, based on the CLSI EP-6A evaluation, the method has been demonstrated to be linear from 1.3 to 26.7 ng/mL, with 5.0 ng/mL up to 66.8 ng/mL, and 25.0 ng/mL in the interval 66.8 and 267.4 ng/mL  $\alpha$ GST. Linearity determined by diluting of human serum clinical sample spiked with known high endogenous  $\alpha$ GST concentration. The FastPack<sup>®</sup> IP  $\alpha$ GST is linear over the range of 0.5 ng/mL (LOQ) to 200 ng/mL. (% Recovery = [Measured value/Expected Value] X 100).

### Linearity Upon Dilution

A serum sample spiked with high endogenous  $\alpha$ GST concentration was diluted 90% to 0.5% with the serum sample and quantified in triplicate determinations with one FastPack<sup>®</sup>  $\alpha$ GST reagent lot on a single FastPack<sup>®</sup> instrument. Measured values of the dilutions compared to their expected values (determined in the same assay) based on the measured value of the Neat sample and respective dilution. (% Recovery = [Measured value/Expected Value] X 100).

%dilution	Expected $\alpha$ GST, ng/mL	Measured $\alpha$ GST, ng/mL	% Recovery
100.0%	N/A	267.4	N/A
90.0%	240.6	229.2	95.2
75.0%	200.5	173.8	86.7
50.0%	133.7	122.1	91.3
25.0%	66.8	54.8	82.0
12.5%	33.4	29.0	86.8
6.3%	16.7	14.3	85.6
3.0%	8.0	6.8	85.0
1.0%	2.7	2.5	92.6
0.5%	1.3	1.2	92.3

In addition, three serum samples with high endogenous  $\alpha$ GST levels were diluted with  $\alpha$ GST Sample Diluent containing phosphate buffered saline solution and Tween 20. The measured values of the dilutions were compared to their expected values.

Serum 1	Expected (ng/mL)	Measured (ng/mL)	% Recovery
Neat	N/A	>200	N/A
1:160	N/A	68.5	N/A
1:320	34.3	32.9	96.6
1:640	17.1	16.4	95.9
1:1280	8.6	7.9	91.9
1:2560	4.3	4.2	97.7


Serum 2	Expected (ng/mL)	Measured (ng/mL)	% Recovery
Neat	N/A	>200	N/A
1:320	N/A	99.9	N/A
1:640	50.0	52.0	104.0
1:1280	25.0	23.9	95.6
1:2560	12.5	12.2	97.6
1:5120	6.25	5.4	86.4

Serum 3	Expected (ng/mL)	Measured (ng/mL)	% Recovery
Neat	N/A	>200	N/A
1:20	N/A	148.9	N/A
1:40	74.4	70.7	95.0
1:80	37.2	35.8	96.2
1:160	18.6	18.2	87.8
1:320	9.3	9.1	97.8

### Method Comparison

Clinical serum samples (n=141) were used to compare the values obtained using the FastPack® IP αGST method and the values obtained using the Argutus Medical HEPKIT® Alpha method. The values were evaluated for agreement using Deming regression analysis, with the associated correlation coefficient.

n	Range of Observation (ng/mL)	Intercept (ng/mL)	Slope	r
141	0.0 – 193.9	1.7	0.85	0.986

### INTERFERING SUBSTANCES

The effect of icteric, lipemic, and hemolyzed serum samples on quantification of αGST was investigated by spiking one serum pool with an αGST concentration of ~60.0 ng/mL with known concentrations of bilirubin, triglycerides, and hemoglobin, and another serum pool with endogenous αGST concentration of ~3.0 ng/mL with the same concentrations of the interfering substances. The values obtained for the serum pools with each interfering substance were compared to the values obtained for the serum pools without the interfering substance and the percentage bias or absolute bias in ng/mL determined. These compounds did not show interference at the levels indicated.

	Interfering Substance		
	Bilirubin (30 mg/dL)	Lipemia (1000 mg/dL)	Hemoglobin (500 mg/dL)
Non-spiked aliquot	3.8 ng/mL	3.0	3.9
Spiked aliquot	3.9 ng/mL	3.2	3.8
Absolute bias, ng/mL	+0.1	+0.2	-0.1

	Interfering Substance		
	Bilirubin (30 mg/dL)	Lipemia (1000 mg/dL)	Hemoglobin (500 mg/dL)
Non-spiked aliquot	66.5 ng/mL	57.7 ng/mL	58.7
Spiked aliquot	66.8 ng/mL	56.2 ng/mL	56.2
% Bias	+0.5%	- 2.6%	- 4.3%

**Cross-reactivity**

No significant cross reactivity (<0.5 ng/mL) was obtained when a FastPack®  $\alpha$ GST sample was spiked with the following cross-reactants: m $\mu$  GST, at 250 ng/mL, and pi GST at 2000 ng/mL.

**Limit of Blank (LOB), Limit of Detection (LOD), and Limit of Quantitation (LOQ)**

The limit of blank (LOB, the highest measurement likely to be observed for a blank sample), limit of detection (LOD, the lowest amount of analyte in a sample that can be detected with type I and II error rates set to 5%), and limit of quantitation (LOQ, the lowest amount of analyte in a sample that can be reliably detected and at which the total error meets the pre-specified requirement for accuracy) were determined according to CLSI EP17-A for the FastPack® IP  $\alpha$ GST assay. In this study, the limit of blank was determined from 20 replicate determinations of the FastPack®  $\alpha$ GST sample on each of three different FastPack® IP instruments. Raw RLUs from the assays were converted to apparent ng/mL based on the calibration curve for each assay. The LOB was determined as the maximum observed value. This value was 0.12 ng/mL  $\alpha$ GST.

The LOD was estimated from 60 replicate determinations of a low control sample per the CLSI EP17-A guideline. In this study, the LOD was found to be 0.21 ng/mL  $\alpha$ GST.

For the LOQ analyses, the prospectively defined goals for accuracy of  $\leq 20\%$  CV. In this study, a control sample precision shown to meet this criterion for CV, thus the LOQ was found to be 0.5 ng /mL  $\alpha$ GST.

## REFERENCES

1. Trull, A.K. *et al.* (1994). Serum  $\alpha$ -Glutathione S-Transferase- A sensitive marker of hepatocellular damage associated with acute liver allograft rejection. *Transplantation*, 58 (12), 1345-1351.
2. Kobayashi, H. *et al.* (2000).  $\alpha$ -Glutathione-S-Transferase as a new sensitive marker of hepatocellular damage in biliary atresia. *Pediatr. Surg. Int.*, 16: 302-305.
3. Mulder, T. P.J. *et al.* (1999). Variability of Glutathione S-Transferase a in human liver and plasma. *Clinical Chemistry* 45 (3), 355-359.
4. Rees, G.W. *et al.* (1995). Evaluation of an enzyme-immunometric assay for serum  $\alpha$ -Gultathione S-Transferase. *Ann. Clin. Biochem.* 32, 575-583.
5. Hughes, V.F. *et al.* (1997). Randomized Trial to evaluate the clinical benefits of serum  $\alpha$ -Glutathione S-Transferase concentration monitoring after liver transplantation. *Transplantation* 64 (10), 1446-1452.
6. Beckett, G.J. and Hayes, J.D. (1993). Glutathione S-Transferases: biomedical applications. *Adv. Clin. Chem.*, 30, 281-389.
7. Platz K.P. *et al.* (1997). Determination of alpha- and Pi-Glutathione-S-Transferase will improve monitoring after liver transplantation. *Transplant Proc.* 1997 Nov; 29(7):2827-9.
8. Schroff, R.J. *et al.* Human anti-mouse immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res*, 45:879 – 885, 1985.
9. Procedures for the collection of diagnostic blood specimens by venipuncture. Approved Standard – Sixth Edition: H3-A6: 27(26), 2007, Clinical and Laboratory Standards Institute (CLSI).
10. Procedures for the handling and processing of blood specimens; Approved guideline – Third Edition, H18-A3; 23 (38), 2004. Clinical and Laboratory Standards Institute (CLSI).

Qualigen, Inc.  
Carlsbad, CA 92011 USA  
Technical Support  
(760) 918-9165  
(877) 709-2169

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## REFERENZEN

1. Trull, A.K. *et al.* (1994). Serum  $\alpha$ -Glutathione S-Transferase- A sensitive marker of hepatocellular damage associated with acute liver allograft rejection. *Transplantation*, 58 (12), 1345-1351.
2. Kobayashi, H. *et al.* (2000).  $\alpha$ -Glutathione-S-Transferase as a new sensitive marker of hepatocellular damage in biliary atresia. *Pediatr. Surg. Int.*, 16: 302-305.
3. Mulder, T. P.J. *et al.* (1999). Variability of Glutathione S-Transferase a in human liver and plasma. *Clinical Chemistry* 45 (3), 355–359.
4. Rees, G.W. *et al.* (1995). Evaluation of an enzyme-immunometric assay for serum  $\alpha$ -Gultathione S-Transferase. *Ann. Clin. Biochem.* 32, 575-583.
5. Hughes, V.F. *et al.* (1997). Randomized Trial to evaluate the clinical benefits of serum  $\alpha$ -Glutathione S-Transferase concentration monitoring after liver transplantation. *Transplantation* 64 (10), 1446-1452.
6. Beckett, G.J. and Hayes, J.D. (1993). Glutathione S-Transferases: biomedical applications. *Adv. Clin. Chem.*, 30, 281–389.
7. Platz K.P. *et al.* (1997). Determination of alpha- and Pi-Glutathione-S-Transferase will improve monitoring after liver transplantation. *Transplant Proc.* 1997 Nov; 29(7):2827-9.
8. Schroff, R.J. *et al.* Human anti-mouse immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res*, 45:879 – 885, 1985.
9. Procedures for the collection of diagnostic blood specimens by venipuncture. Approved Standard – Sixth Edition: H3-A6: 27(26), 2007, Clinical and Laboratory Standards Institute (CLSI).
10. Procedures for the handling and processing of blood specimens; Approved guideline – Third Edition, H18-A3; 23 (38), 2004. Clinical and Laboratory Standards Institute (CLSI).

Qualigen, Inc.  
Carlsbad, CA 92011 USA  
Technischer Kundendienst  
+1 (760) 918-9165  
+1 (877) 709-2169

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