The FastPack® IP Free PSA Assay should be used only with the FastPack® IP and FastPack® Total PSA Assay to calculate the ratio of free PSA to total PSA (percent free PSA). Results contained in this insert apply only to percent free PSA as measured by the FastPack® IP Free PSA and Total PSA Assays.

The concentration of free PSA and total PSA in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the free and total PSA assay method used. Values obtained using different assay methods cannot be used interchangeably.

Free PSA concentrations are dependent on the standard used to calibrate the assay. Free PSA concentrations based on calibration to the WHO 96/670 Reference Preparation will differ significantly from Free PSA concentrations based on calibration to the original Hybritech assay. The concentrations are not inter-changeable. If the calibration is changed, accepted laboratory practice is to establish a new baseline for patient monitoring.¹ For more information on WHO and Hybritech configuration reference your instrument procedure manual.

CAUTION: United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.

INTENDED USE

The FastPack® IP Free PSA Immunoassay is a paramagnetic particle immunoassay for the in vitro quantitative determination of free prostate-specific antigen (free PSA) in human serum. The FastPack® IP Free PSA Immunoassay is designed for use in conjunction with the FastPack® IP and FastPack® Total PSA Immunoassay and the FastPack® IP System for calculating the ratio of free/total PSA, expressed as a percentage (percent free PSA).

The percent free PSA, as measured by the FastPack® IP System, is intended for use as an aid to differentiate between prostate cancer and benign prostatic conditions, in male patients over 50 with total PSA values between 4 and 10 ng/ml and digital rectal examinations that are not suspicious for cancer. Prostatic biopsy is required for diagnosis of cancer.

SUMMARY

Prostate-specific antigen (PSA) is a glycoprotein (molecular weight 30,000-34,000 Daltons), which exhibits a high degree of homology with serine proteases of the kallikrein family. It has the function of serine protease with chymotrypsin-like activity.² The proteolytic activity of PSA in blood is inhibited by the irreversible formation of complexes with protease inhibitors such as α -1- antichymotrypsin (ACT), α -2-macroglobulin and other acute phase proteins.³ PSA occurs in three forms in blood: two immunodetectable forms are PSA complexed with the serine protease inhibitor α -1-antichymotrypsin; the second form is free or uncomplexed PSA.^{4, 5, 6} The third form that is undetectable in conventional immunoassay is PSA complexed with α -2 macroglobulin. This is due to engulfment and masking of PSA epitopes by the α -2 macroglobulin molecule.

Elevated levels of PSA in serum are generally indicative of a pathological condition of the prostate, for example, carcinoma, benign prostatic hyperplasia (BPH) or prostatitis.^{7,8,9} Prostate Specific Antigen is also present in paraurethral and anal glands, as well as breast tissue. An inflammation, trauma or stimulation of the prostate (e.g. following digital rectal examination, biopsy and colonoscopy, etc.) can result in PSA elevations of varying duration and magnitude. It has been demonstrated that the proportion of PSA-ACT in serum is higher in prostate cancer patients than in normal subjects and/or in patients with noncancerous prostatic diseases such as BPH or prostatitis.¹⁰

Determination of the free PSA level in a patient's sample and calculation of the percentage of free to Total PSA enhances the ability to differentiate between prostate cancer and prostatic diseases such as BPH. Measurement of free PSA thus reduces unnecessary biopsies and results in a reduction in cost.

TEST PRINCIPAL

The FastPack® IP Free PSA Immunoassay is a chemiluminescence assay based on the "sandwich" principle.

- Primary incubation: Specimen, control or calibrator [100 μ L] and antibody solution (mixture of a biotinylated monoclonal PSA-specific antibody and a monoclonal free PSA-specific antibody labeled with alkaline phosphatase) [100 μ L] react to form a sandwich complex.
- Secondary incubation: Streptavidin-coated paramagnetic particle solution (solid-phase) is added to the reaction mixture. During this incubation, the sandwich complex is bound to the solid-phase via the interaction of biotin and streptavidin.
- Removal of unbound materials: The paramagnetic particles are washed with wash buffer [0.2 mL/wash] to remove unbound materials.

- Substrate addition and detection: Chemiluminogenic substrate [140 µL] is added to the solid-phase bound complex and results in “glow” chemiluminescence, which is measured using the FastPack® IP Analyzer.
- The amount of bound labeled-antibody is directly proportional to the concentration of free PSA in the sample.

REAGENTS - Content and Concentration

FastPack IP® Free PSA Immunoassay Kit – Cat. No. 25000045

Each FastPack® IP carton contains:

- 30 FastPack® IPs

Each FastPack® IP Contains:

- Paramagnetic Particles, 150 µL
Streptavidin-coated paramagnetic particles in buffer containing 0.1% sodium azide as a preservative.
- Free PSA Antibody Solution, 100 µL
Antibody solution containing mouse monoclonal antibody coupled to biotin and mouse monoclonal antibody labeled with alkaline phosphatase in a protein matrix containing 0.1% sodium azide as preservative.
- Wash Buffer, 2.0 mL
Tris buffer containing surfactants.
- Substrate, 140 µL
ImmuGlow™: Indoxyl-3-phosphate and lucigenin in buffer containing preservatives.

Materials required but not provided

- FastPack® IP System
- FastPack® Free PSA Calibrator Kit - Cat. No. 25000008
- FastPack® PSA Control Kit - Cat. No. 25000003 (International)

WARNING AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Do not eat, drink or smoke in designated work areas.
- Wash hands thoroughly after handling specimen.
- HAMA Interference: some individuals have antibodies to mouse protein (HAMA), which can cause interference in immunoassays that employ antibodies derived from mice. In particular, it has been reported that serum samples from patients who have undergone therapy or diagnostic procedures that include infusion of mouse monoclonal antibody may produce erroneous results in such assays.
- FastPack® IP reagents are stable until the expiration date on the label when stored and handled as directed. Do not use FastPack® IP reagents beyond expiration date.
- Discard used FastPacks into a biohazard container.
- The components containing Sodium Azide are classified per applicable European Economic Community (EEC) Directives as: Harmful (Xn). The following are appropriate Risk (R) and Safety (S) phrases:

R22 Harmful if swallowed.
 R32 Contact with acids liberates very toxic gas.
 S2 Keep out of reach of children.
 S13 Keep away from food, drink and animal feeding stuffs.
 S36 Wear suitable protective clothing.
 S46 If swallowed, seek medical advice immediately and show this container or label.

STORAGE INSTRUCTIONS

Store at 2 - 8 °C. Protect from light.

SPECIMEN COLLECTION/PREPARATION

1. Serum is required for the FastPack® IP Free PSA Immunoassay. Plasma samples should not be used.
2. The National Committee for Clinical Laboratory Standards (NCCLS) provides the following recommendations for handling, processing, and storing blood.^{11,12}
 - Collect all blood samples observing routine precautions for venipuncture.
3. It is not required that patients fast prior to blood collection.
4. For serum samples:
 - Ensure that complete clot formation has occurred before centrifugation. This takes approximately 30 minutes. Some samples may exhibit increased clotting time, especially those from patients receiving anticoagulant or thrombolytic therapy.
 - Serum should be centrifuged and separated from the clot within 3 hours from time of collection.
 - Remove serum from the cells prior to storage at 2-8 °C.
 - If not tested within 24 hours, the sample should be frozen at -20 °C or colder.
5. Do not freeze samples (-20 °C) for more than two months.
6. Frozen samples should be thawed completely and mixed by gentle inversion prior to use.
7. Samples should be free of Fibrin, Red Blood Cells, or other particulate material for optimal results. Serum samples showing turbidity and/or particulate matter should be centrifuged prior to use.
8. Ensure the samples are free of bubbles.

9. Samples collected up to 2 hours after DRE show no significant increases in PSA.¹³
10. Human samples should be handled in accordance with the OSHA standard on Bloodborne Pathogens.²¹

ASSAY PROCEDURE

See the FastPack® IP System Procedure Manual for information on operating the FastPack® IP System.

INSTRUMENTATION

FastPack® IP System

DETAILS OF CALIBRATION

During the FastPack® IP production process, Qualigen generates a master standard curve and places this information in the barcode of each FastPack® IP label, where it can be read by the FastPack® IP Analyzer during the testing sequence. The FastPack® IP Analyzer must be calibrated by the user to ensure that it is properly adjusted for the particular lot of FastPacks that are being used. Separate calibrations must be run for each type of test, i.e. Free PSA, Total PSA or Testosterone. The frequency of calibration varies for each test type. For the FastPack® IP Free PSA Immunoassay, the FastPack® IP Analyzer must be calibrated once every 14 days or whenever a new lot of Free PSA FastPacks are to be used.

Whenever the user performs an initial calibration for a particular lot of FastPacks or uses a new lot of calibrator, 2 FastPacks must be run for calibration (duplicates). Whenever recalibration is performed with the same lot of FastPacks and calibrator, 2 FastPacks must be run for calibration. See FastPack® IP System Procedure Manual for "Running a Calibration".

Use FastPack® Free PSA Calibrator Kit – Cat. No. 25000008

RESULTS

The FastPack® IP Analyzer uses the information from the barcode to construct a lookup table of x,y values that represent the standard curve and estimates the concentration of unknown samples by linear interpolation.

QUALITY CONTROL

Quality control materials simulate real specimens and are essential for monitoring the system performance of assays. Good Laboratory Practices (GLP) include the use of control specimens to ensure that all reagents and protocols are performing properly. See FastPack® IP System Procedure Manual for "Running Controls".

Controls available: FastPack® PSA Control Kit – Cat. No. 25000003 (International)

LIMITATION OF PROCEDURE

- The FastPack® IP Free PSA is not intended to be used to diagnose prostate cancer.
- Samples can be accurately measured within the reportable range of the analytical sensitivity and the highest calibrator, 20 ng/mL (16.0 ng/mL WHO).
- Samples > 20 ng/mL (16.0 ng/mL WHO) should be run using another method. Dilution of out of range results is not recommended.
- FastPack® IP Free PSA Immunoassay is tested up to 1000 ng/mL and showed no hook effect.
- Samples from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such samples may show either falsely elevated or depressed values when tested with assay kits that employ mouse monoclonal antibodies.^{14,15}
- It is known that in rare cases PSA isoforms do exist which may be measured differently by different PSA tests. Findings of this kind have occasionally been reported for PSA tests from various manufacturers.^{16,17,18}
- Heterophilic antibodies in serum have the potential to cause interference in immunoassay systems.^{19,20} Infrequently, PSA levels may appear elevated due to heterophilic antibodies present in the patient's serum or to nonspecific protein binding. If the PSA level is inconsistent with clinical evidence, additional PSA testing is suggested to confirm the result.
- Prostatic massage, ultrasonography and needle biopsy may cause clinically significant elevation in PSA levels. Hormonal therapy may affect PSA expression, therefore a low PSA level after treatment that includes hormonal therapy may not adequately reflect the presence of residual or recurrent disease.
- The results from the FastPack® IP Free PSA Immunoassay should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

CLINICAL SIGNIFICANCE

Prostate cancer is more likely to produce bound PSA. Therefore, the ratio of free to total PSA is used to determine how much PSA exists as bound PSA. The relative risk of prostate cancer is as low as 8% if the free PSA ratio is over 25% and as high as 56% if the ratio is <10%. In Caucasian men with total PSA between 4 & 10 ng/mL, 95% of those found to have cancer by biopsy have a free to total ratio of <24% and 90% have a ratio of <22%. Based on these results, some urologists recommend prostate biopsy for men with ratios <25%. Men with higher ratios can be followed with annual examinations. Using a single decision point of 25% free PSA as an indication for biopsy will detect 95% of cancers and avoid 25% of unnecessary biopsies. Free PSA ratio can also be used to counsel patients about their individual cancer risk. Patients can then decide to have a biopsy or wait and reevaluate PSA in 6 months.²¹ Note: referenced study was not using the FP system and is only presented as a basis for scientific validity.

A total of 101 patient samples with the 100µL tPSA assay and 100 samples with the 25µL tPSA assay recovered between the 4-10 ng/mL. The table below represents a distribution of results that fell within each of the Free to Total PSA ratios at a certain age using the FP system and do not represent cancer probability.

Total PSA 100 µL

PSA (ng/mL)	Free to Total ratio (%)	Age (years)		
		50-59	60-69	70+
4 to 10	<10	36.4% (8)	11.1% (5)	20.6% (7)
4 to 10	10-25	54.5% (12)	68.9% (31)	64.7% (22)
4 to 10	>25	9.1% (2)	20% (9)	14.7% (5)

Total PSA 25 µL

PSA (ng/mL)	Free to Total ratio (%)	Age (years)		
		50-59	60-69	70+
4 to 10	<10	31.8% (7)	17.4% (8)	25% (8)
4 to 10	10-25	63.6% (14)	65.2% (30)	65.6% (21)
4 to 10	>25	4.5% (1)	17.4% (8)	9.4% (3)

Total PSA 100 µL

WHO PSA (ng/mL)	Free to Total ratio (%)	Age (years)		
		50-59	60-69	70+
3.2 to 8	<10	36.4% (8)	11.1% (5)	20.6% (7)
3.2 to 8	10-25	54.5% (12)	68.9% (31)	64.7% (22)
3.2 to 8	>25	9.1% (2)	20% (9)	14.7% (5)

Total PSA 25 µL

WHO PSA (ng/mL)	Free to Total ratio (%)	Age (years)		
		50-59	60-69	70+
3.2 to 8	<10	31.8% (7)	17.4% (8)	25% (8)
3.2 to 8	10-25	63.6% (14)	65.2% (30)	65.6% (21)
3.2 to 8	>25	4.5% (1)	17.4% (8)	9.4% (3)

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The reproducibility of the Free PSA assay was measured by assaying five different levels (n=80, for each level) twenty non-consecutive days using two analyzers and two lots of reagents. The coefficient of variation (%CV) between analyzers, between runs was calculated using analysis of variance.

Lot 1

Sample	Average, ng/mL		Within Run	Between Run	Between Day	Total Precision	%TE	%T'Ea
	Hybritech	WHO	%CV	%CV	%CV	%CV		
1	0.34	0.27	10.416	2.812	6.019	12.354	27.454	30
2	1.84	1.47	7.065	0.000	5.761	9.116	28.344	30
3	5.55	4.44	6.554	2.725	6.618	9.705	17.550	30
4	10.86	8.69	6.770	0.000	7.026	9.757	21.756	30
5	16.60	13.28	5.434	1.852	8.890	10.582	23.599	30

Lot 2

Sample	Average, ng/mL		Within Run	Between Run	Between Day	Total Precision	%TE	%T'Ea
	Hybritech	WHO	%CV	%CV	%CV	%CV		
1	0.44	0.35	5.794	4.492	0.534	7.351	19.438	30
2	2.28	1.82	5.241	0.000	2.078	5.638	10.917	30
3	6.16	4.93	4.527	0.000	1.361	4.727	12.199	30
4	11.19	8.96	4.519	2.102	0.516	5.011	11.174	30
5	16.17	12.93	5.520	0.000	2.482	6.053	11.688	30

Method Comparison

152 serum samples were used to compare FastPack® Free PSA against predicate Beckman Access Hybritech free PSA. Evaluation was performed using Passing-Bablok regression analysis, with a slope of 1.09, intercept of 0.006 ng/mL (Hybritech) and r-value of 0.986.

N	Range of Values (ng/mL)	Intercept (ng/mL)	Slope	R
152	0.16 to 17.2	0.006	1.09	0.986

N	Range of Values, WHO (ng/mL)	Intercept (ng/mL)	Slope	R
152	0.13 to 13.8	0.005	1.09	0.986

INTERFERING SUBSTANCES

Interfering substances were added to serum pools or controls containing known amounts of free PSA. The value obtained for the sample with each interfering substance was compared to the value obtained for the sample without the interfering substance. The following compounds did not show interference at the levels indicated.

Test Compound	Test Concentration	Chemotherapeutic Agents	Concentration
d-Biotin	500 ng/mL	Cyclophosphamide	700 µg/mL
Bilirubin	49 mg/dL	Diethylstilbestrol	2 µg/mL
Hemoglobin	600 mg/dL	Doxorubicin HCl	16 µg/mL
Human IgG	1000 mg/dL	Methotrexate	8 µg/mL
Prostatic Acid Phosphatase (PAP)	1 µg/dL	Lupron	100 µg/mL
Triglycerides	1000 mg/dL	Hytrin	40 µg/mL
Total Protein	5 gm/dL	Proscar	100 µg/mL
		Ketoconazole	200 µg/mL

* Interferences was found with d-Biotin above 500 ng/mL.

Analytical Sensitivity:**Limit of Blank (LOB), Limit of Detection (LOD), and Limit of Quantitation (LOQ)**

The limit of blank (LOB, the highest measurement likely to be observed for a blank sample), limit of detection (LOD, the lowest amount of analyte in a sample that can be detected with type I and II error rates set to 5%), and limit of quantitation (LOQ, the lowest amount of analyte in a sample that can be reliably detected) were determined according to CLSI EP17-A2. In this study the limit of blank was determined from 240 replicate determinations of a blank sample tested on six different FastPack® analyzers using three reagent lots. Raw RLUs from the assays were converted to apparent ng/mL based on the calibration curve for each assay. The LOB was determined as the 228th rank of the sorted distribution of values. This value was 0.01 ng/mL FPSA.

The LOD was estimated from 180 replicate determinations of four low level samples. Per the CLSI EP17-A2 guideline, the parametric LOD calculation was utilized and yielded 0.03 ng/mL FPSA.

The LOQ was determined as the lowest sample which provided <20% CV. The LOQ was set to 0.05 ng/mL.

Analytical Specificity

Analytical Specificity was measured by spiking PSA/ACT complex into Calibrator Matrix at levels of 1000 ng/mL, 500 ng/mL, and 100 ng/mL. Less than 1% cross-reactivity was observed in each case.

REFERENCES

- ¹ Lilja H, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines (LMPG): Practice Guidelines and Recommendations for Use of Tumor Markers in the Clinic, Prostate Cancer (Section B), Draft 2006. National Academy of Clinical Biochemistry.
- ¹ Henttu P, Vihko P. Prostate-specific antigen and human glandular kallikrein: Two kallikreins of the human prostate. *Ann Med* 1994;26:157-164.
- ² Tewari PC, Bluestein BI. Multiple forms of prostate specific antigen and the influences of immunoassay design on their measurement in patient serum. *J Clin Ligand Assay*, 18 1995;3:186-196.
- ³ Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and characterization of different molecular forms of prostatespecific antigen in human seminal fluid. *Clin Chem*. 1995;41(11):1567-1573.
- ⁴ Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. *Clin Lab Invest* 1995;55 Suppl221:32-34.
- ⁵ Scher HI, Kelly WK. Flutamide withdrawal syndrome: Its impact on clinical trials in hormone-refractory prostate cancer. *J Clin Oncol* 1993;11(8):1566-1572.
- ⁶ Partin AW, et al. Prostate specific antigen in the staging of localized prostate cancer: Influence of tumor differentiation, tumor volume and benign hyperplasia. *J Urol* 1990; 143:747-752.
- ⁷ Semjonow A, Brandt B, Oberpenning F, Hertle L. Discrepancies in assays impair the interpretation of prostate-specific antigen. *Urology* 1995;34:303-315.
- ⁸ Armbruster DA. Prostate-Specific Antigen: Biochemistry, Analytical, and Clinical Application [Review]. *Clin Chem* 1993; 39:181.
- ⁹ Christenson A, Bjork T Nilsson O, Dakleen U, Matikainen MT, Cockett ATK, Abrahamsson PA to α -1- antichymotrypsin and an indicator of Prostate Cancer. *J. Urol* 150: 100-105, 1993
- ¹⁰ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition; H3-A5; 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS).
- ¹¹ Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- ¹² Chybowski FM, Berstralh EJ, Oesterling JE. The effect of digital rectal examination on the serum prostate specific antigen concentration: Results of a randomized study. *J Urol* 1992;148:83-86.
- ¹³ Primus FJ, Kelly EA, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem*. 1988;34(2):261-264.
- ¹⁴ Schroff RJ, Foon KA, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-885.
- ¹⁵ Van Duijnhoven HLP, Pérquériaux NCV, van Zon JPHM, Blankenstein MA. Large discrepancy between prostate-specific antigen results from different assays during longitudinal follow-up of a prostate cancer patient. *Clin Chem* 1996;42(4):637-641.
- ¹⁶ Wians FH. The "correct" PSA concentration. *Clin Chem*. 1996;42(11):1882-1885.
- ¹⁷ Cohen RJ, Haffejee Z, Steele GS, Naylar SJ. Advanced prostate cancer with normal serum prostate-specific antigen levels. *Arch Pathol Lab Med* 1994;118:1123-1126.
- ¹⁸ Nahm MH, Hoffman JW. Heteroantibody: Phantom of the immunoassay. *Clin Chem* 1990;36(6):829.
- ¹⁹ Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- ²⁰ US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. *Federal Register* 1991; 56(235): 64175-82.
- ²¹ Catalona WJ, Partin AW, Slawin KM, et al. Use of the Percentage of Free Prostate-Specific Antigen to Enhance Differentiation of Prostate Cancer From Benign Prostatic Disease: A Prospective Multicenter Clinical Trial. *JAMA*. 1998;279(19):1542–1547. doi:10.1001/jama.279.19.1542

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REFERENZEN

- ¹ Lilja H, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines (LMPG): Practice Guidelines and Recommendations for Use of Tumor Markers in the Clinic, Prostate Cancer (Section B), Draft 2006. National Academy of Clinical Biochemistry.
- ² Henttu P, Vihko P. Prostate-specific antigen and human glandular kallikrein: Two kallikreins of the human prostate. *Ann Med* 1994;26:157-164.
- ³ Tewari PC, Bluestein BI. Multiple forms of prostate specific antigen and the influences of immunoassay design on their measurement in patient serum. *J Clin Ligand Assay*, 18 1995;3:186-196.
- ⁴ Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and characterization of different molecular forms of prostatespecific antigen in human seminal fluid. *Clin Chem*. 1995;41(11):1567-1573.
- ⁵ Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. *Clin Lab Invest* 1995;55 Suppl221:32-34.
- ⁶ Scher HI, Kelly WK. Flutamide withdrawal syndrome: Its impact on clinical trials in hormone-refractory prostate cancer. *J Clin Oncol* 1993;11(8):1566-1572.
- ⁷ Partin AW, et al. Prostate specific antigen in the staging of localized prostate cancer: Influence of tumor differentiation, tumor volume and benign hyperplasia. *J Urol* 1990; 143:747-752.
- ⁸ Semjonow A, Brandt B, Oberpenning F, Hertle L. Discrepancies in assays impair the interpretation of prostate-specific antigen. *Urology* 1995;34:303-315.
- ⁹ Armbruster DA. Prostate-Specific Antigen: Biochemistry, Analytical, and Clinical Application [Review]. *Clin Chem* 1993; 39:181.
- ¹⁰ Christenson A, Bjork T Nilsson O, Dakleen U, Matikainen MT, Cockett ATK, Abrahamsson PA to α -1- antichymotrypsin and an indicator of Prostate Cancer. *J. Urol* 150: 100-105, 1993
- ¹¹ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition; H3-A5; 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS).
- ¹² Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- ¹³ Chybowski FM, Berstralh EJ, Oesterling JE. The effect of digital rectal examination on the serum prostate specific antigen concentration: Results of a randomized study. *J Urol* 1992;148:83-86.
- ¹⁴ Primus FJ, Kelly EA, et al. „Sandwich“-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem*. 1988;34(2):261-264.
- ¹⁵ Schroff RJ, Foon KA, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-885.
- ¹⁶ Van Duijnhoven HLP, Pérquériaux NCV, van Zon JPHM, Blankenstein MA. Large discrepancy between prostate-specific antigen results from different assays during longitudinal follow-up of a prostate cancer patient. *Clin Chem* 1996;42(4):637-641.
- ¹⁷ Wians FH. The „correct“ PSA concentration. *Clin Chem*. 1996;42(11):1882-1885.
- ¹⁸ Cohen RJ, Haffejee Z, Steele GS, Naylar SJ. Advanced prostate cancer with normal serum prostate-specific antigen levels. *Arch Pathol Lab Med* 1994;118:1123-1126.
- ¹⁹ Nahm MH, Hoffman JW. Heteroantibody: Phantom of the immunoassay. *Clin Chem* 1990;36(6):829.
- ²⁰ Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- ²¹ US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. *Federal Register* 1991; 56(235): 64175-82.
- ²¹ Catalona WJ, Partin AW, Slawin KM, et al. Use of the Percentage of Free Prostate-Specific Antigen to Enhance Differentiation of Prostate Cancer From Benign Prostatic Disease: A Prospective Multicenter Clinical Trial. *JAMA*. 1998;279(19):1542–1547. doi:10.1001/jama.279.19.1542

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RÉFÉRENCES

- ¹ Lilja H, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines (LMPG): Practice Guidelines and Recommendations for Use of Tumor Markers in the Clinic, Prostate Cancer (Section B), Draft 2006. National Academy of Clinical Biochemistry.
- ² Henttu P, Vihko P. Prostate-specific antigen and human glandular kallikrein: Two kallikreins of the human prostate. *Ann Med* 1994;26:157-164.
- ³ Tewari PC, Bluestein BI. Multiple forms of prostate specific antigen and the influences of immunoassay design on their measurement in patient serum. *J Clin Ligand Assay*, 18 1995;3:186-196.
- ⁴ Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and characterization of different molecular forms of prostatespecific antigen in human seminal fluid. *Clin Chem*. 1995;41(11):1567-1573.
- ⁵ Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. *Clin Lab Invest* 1995;55 Suppl221:32-34.
- ⁶ Scher HI, Kelly WK. Flutamide withdrawal syndrome: Its impact on clinical trials in hormone-refractory prostate cancer. *J Clin Oncol* 1993;11(8):1566-1572.
- ⁷ Partin AW, et al. Prostate specific antigen in the staging of localized prostate cancer: Influence of tumor differentiation, tumor volume and benign hyperplasia. *J Urol* 1990; 143:747-752.
- ⁸ Semjonow A, Brandt B, Oberpenning F, Hertle L. Discrepancies in assays impair the interpretation of prostate-specific antigen. *Urology* 1995;34:303-315.
- ⁹ Armbruster DA. Prostate-Specific Antigen: Biochemistry, Analytical, and Clinical Application [Review]. *Clin Chem* 1993; 39:181.
- ¹⁰ Christenson A, Bjork T Nilsson O, Dakleen U, Matikainen MT, Cockett ATK, Abrahamsson PA to α -1- antichymotrypsin and an indicator of Prostate Cancer. *J. Urol* 150: 100-105, 1993
- ¹¹ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition; H3-A5; 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS).
- ¹² Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- ¹³ Chybowski FM, Berstralh EJ, Oesterling JE. The effect of digital rectal examination on the serum prostate specific antigen concentration: Results of a randomized study. *J Urol* 1992;148:83-86.
- ¹⁴ Primus FJ, Kelly EA, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem*. 1988;34(2):261-264.
- ¹⁵ Schroff RJ, Foon KA, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-885.
- ¹⁶ Van Duijnhoven HLP, Pérquériaux NCV, van Zon JPHM, Blankenstein MA. Large discrepancy between prostate-specific antigen results from different assays during longitudinal follow-up of a prostate cancer patient. *Clin Chem* 1996;42(4):637-641.
- ¹⁷ Wians FH. The "correct" PSA concentration. *Clin Chem*. 1996;42(11):1882-1885.
- ¹⁸ Cohen RJ, Haffejee Z, Steele GS, Naylar SJ. Advanced prostate cancer with normal serum prostate-specific antigen levels. *Arch Pathol Lab Med* 1994;118:1123-1126.
- ¹⁹ Nahm MH, Hoffman JW. Heteroantibody: Phantom of the immunoassay. *Clin Chem* 1990;36(6):829.
- ²⁰ Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- ²¹ US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. *Federal Register* 1991; 56(235): 64175-82.
- ²¹ Catalona WJ, Partin AW, Slawin KM, et al. Use of the Percentage of Free Prostate-Specific Antigen to Enhance Differentiation of Prostate Cancer From Benign Prostatic Disease: A Prospective Multicenter Clinical Trial. *JAMA*. 1998;279(19):1542-1547. doi:10.1001/jama.279.19.1542

Qualigen, Inc.
Carlsbad, CA 92011 États-Unis
Assistance technique:
+1 (760) 918-9165
+1 (877) 709-2169

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REFERENCIAS

- ¹ Lilja H, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines (LMPG): Practice Guidelines and Recommendations for Use of Tumor Markers in the Clinic, Prostate Cancer (Section B), Draft 2006. National Academy of Clinical Biochemistry.
- ² Henttu P, Vihko P. Prostate-specific antigen and human glandular kallikrein: Two kallikreins of the human prostate. *Ann Med* 1994;26:157-164.
- ³ Tewari PC, Bluestein BI. Multiple forms of prostate specific antigen and the influences of immunoassay design on their measurement in patient serum. *J Clin Ligand Assay*, 18 1995;3:186-196.
- ⁴ Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and characterization of different molecular forms of prostatespecific antigen in human seminal fluid. *Clin Chem*. 1995;41(11):1567-1573.
- ⁵ Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. *Clin Lab Invest* 1995;55 Suppl221:32-34.
- ⁶ Scher HI, Kelly WK. Flutamide withdrawal syndrome: Its impact on clinical trials in hormone-refractory prostate cancer. *J Clin Oncol* 1993;11(8):1566-1572.
- ⁷ Partin AW, et al. Prostate-specific antigen in the staging of localized prostate cancer: Influence of tumor differentiation, tumor volume and benign hyperplasia. *J Urol* 1990; 143:747-752.
- ⁸ Semjonow A, Brandt B, Oberpenning F, Hertle L. Discrepancies in assays impair the interpretation of prostate-specific antigen. *Urology* 1995;34:303-315.
- ⁹ Armbruster DA. Prostate-Specific Antigen: Biochemistry, Analytical, and Clinical Application [Review]. *Clin Chem* 1993; 39:181.
- ¹⁰ Christenson A, Bjork T Nilsson O, Dakleen U, Matikainen MT, Cockett ATK, Abrahamsson PA to α -1- antichymotrypsin and an indicator of Prostate Cancer. *J. Urol* 150: 100-105, 1993
- ¹¹ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition; H3-A5; 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS).
- ¹² Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- ¹³ Chybowski FM, Berstralh EJ, Oesterling JE. The effect of digital rectal examination on the serum prostate specific antigen concentration: Results of a randomized study. *J Urol* 1992;148:83-86.
- ¹⁴ Primus FJ, Kelly EA, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem*. 1988;34(2):261-264.
- ¹⁵ Schroff RJ, Foon KA, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-885.
- ¹⁶ Van Duijnhoven HLP, Pérquériaux NCV, van Zon JPHM, Blankenstein MA. Large discrepancy between prostate-specific antigen results from different assays during longitudinal follow-up of a prostate cancer patient. *Clin Chem* 1996;42(4):637-641.
- ¹⁷ Wians FH. The "correct" PSA concentration. *Clin Chem*. 1996;42(11):1882-1885.
- ¹⁸ Cohen RJ, Haffejee Z, Steele GS, Naylar SJ. Advanced prostate cancer with normal serum prostate-specific antigen levels. *Arch Pathol Lab Med* 1994;118:1123-1126.
- ¹⁹ Nahm MH, Hoffman JW. Heteroantibody: Phantom of the immunoassay. *Clin Chem* 1990;36(6):829.
- ²⁰ Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- ²¹ US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. *Federal Register* 1991; 56(235): 64175-82.
- ²¹ Catalona WJ, Partin AW, Slawin KM, et al. Use of the Percentage of Free Prostate-Specific Antigen to Enhance Differentiation of Prostate Cancer From Benign Prostatic Disease: A Prospective Multicenter Clinical Trial. *JAMA*. 1998;279(19):1542–1547. doi:10.1001/jama.279.19.1542

Qualigen, Inc.
Carlsbad, CA 92011 EE.
UU.
Soporte técnico:
+1 (760) 918-9165
+1 (877) 709-2169

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BIBLIOGRAFIA

- 1 Lilja H, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines (LMPG): Practice Guidelines and Recommendations for Use of Tumor Markers in the Clinic, Prostate Cancer (Section B), Draft 2006. National Academy of Clinical Biochemistry.
- 2 Henttu P, Vihko P. Prostate-specific antigen and human glandular kallikrein: Two kallikreins of the human prostate. *Ann Med* 1994;26:157-164.
- 3 Tewari PC, Bluestein BI. Multiple forms of prostate specific antigen and the influences of immunoassay design on their measurement in patient serum. *J Clin Ligand Assay*, 18 1995;3:186-196.
- 4 Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and characterization of different molecular forms of prostatespecific antigen in human seminal fluid. *Clin Chem*. 1995;41(11):1567-1573.
- 5 Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. *Clin Lab Invest* 1995;55 Suppl221:32-34.
- 6 Scher HI, Kelly WK. Flutamide withdrawal syndrome: Its impact on clinical trials in hormone-refractory prostate cancer. *J Clin Oncol* 1993;11(8):1566-1572.
- 7 Partin AW, et al. Prostate specific antigen in the staging of localized prostate cancer: Influence of tumor differentiation, tumor volume and benign hyperplasia. *J Urol* 1990; 143:747-752.
- 8 Semjonow A, Brandt B, Oberpenning F, Hertle L. Discrepancies in assays impair the interpretation of prostate-specific antigen. *Urology* 1995;34:303-315.
- 9 Armbruster DA. Prostate-Specific Antigen: Biochemistry, Analytical, and Clinical Application [Review]. *Clin Chem* 1993; 39:181.
- 10 Christenson A, Bjork T Nilsson O, Dakleen U, Matikainen MT, Cockett ATK, Abrahamsson PA to α -1- antichymotrypsin and an indicator of Prostate Cancer. *J. Urol* 150: 100-105, 1993
- 11 Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition; H3-A5; 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS).
- 12 Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- 13 Chybowski FM, Berstralh EJ, Oesterling JE. The effect of digital rectal examination on the serum prostate specific antigen concentration: Results of a randomized study. *J Urol* 1992;148:83-86.
- 14 Primus FJ, Kelly EA, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem*. 1988;34(2):261-264.
- 15 Schroff RJ, Foon KA, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-885.
- 16 Van Duijnhoven HLP, Pérquériaux NCV, van Zon JPHM, Blankenstein MA. Large discrepancy between prostate-specific antigen results from different assays during longitudinal follow-up of a prostate cancer patient. *Clin Chem* 1996;42(4):637-641.
- 17 Wians FH. The "correct" PSA concentration. *Clin Chem*. 1996;42(11):1882-1885.
- 18 Cohen RJ, Haffejee Z, Steele GS, Naylar SJ. Advanced prostate cancer with normal serum prostate-specific antigen levels. *Arch Pathol Lab Med* 1994;118:1123-1126.
- 19 Nahm MH, Hoffman JW. Heteroantibody: Phantom Pf the immunoassay. *Clin Chem* 1990;36(6):829.
- 20 Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- 21 US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. *Federal Register* 1991; 56(235): 64175-82.
- 21 Catalona WJ, Partin AW, Slawin KM, et al. Use of the Percentage of Free Prostate-Specific Antigen to Enhance Differentiation of Prostate Cancer From Benign Prostatic Disease: A Prospective Multicenter Clinical Trial. *JAMA*. 1998;279(19):1542-1547. doi:10.1001/jama.279.19.1542

Qualigen, Inc.
Carlsbad, CA 92011 USA
Assistenza tecnica:
+1 (760) 918-9165
+1 (877) 709-2169

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© 2000 Qualigen, Inc. Alle Rechte vorbehalten. Qualigen und FastPack sind Warenzeichen oder eingetragene Warenzeichen von Qualigen, Inc. Alle sonstigen Warenzeichen sind Eigentum der jeweiligen Besitzer.