The concentration of Thyroid Stimulating Hormone (TSH) in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the TSH assay method used. Values obtained with different assay methods should not be used interchangeably.

**CAUTION:** United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.

### INTENDED USE

The FastPack® IPIPTSH Immunoassay is a paramagnetic particle immunoassay for the *in-vitro* quantitative determination of TSH in human serum and plasma. The measurements of thyroid stimulating hormone (TSH) produced by the anterior pituitary are used in the diagnosis of thyroid or pituitary disorders. The FastPack® IPTSH Immunoassay is designed for use with the FastPack® IP System and the FastPack® System.

### SUMMARY

Thyroid hormones (T4 and T3) have many actions and are indispensable for growth, development, and sexual maturation. Actions include stimulation of heart rate and heart contraction, stimulation of protein synthesis and carbohydrate metabolism, increase in the synthesis and degradation of cholesterol and triglycerides, increase in vitamin requirements, and enhancement of sensitivity of  $\beta$ -adrenergic receptors to catecholamines<sup>1,2</sup>. Therefore, thyroid dysfunction can have a profound global effect.

A tightly coordinated feedback relationship exists among the thyroid gland, hypothalamus, and pituitary gland<sup>3</sup>. The systems of these glands are closely interrelated and integrated, the result being maintenance of thyroid hormone levels in blood. Thyrotropin-releasing hormone (TRH), a tripeptide, is produced in the hypothalamus. TRH acts on the pituitary thyrotropes to cause synthesis and release of TSH<sup>7</sup>. TSH in turn controls the synthesis and release of thyroid hormones. A rise in thyroid hormone levels inhibits the pituitary response to TRH. A fall in thyroid hormone levels causes an increase in TRH and TSH secretion. A failure at any regulatory level of production of thyroid hormones (T4 and T3) will result in either underproduction (hypothyroidism) or overproduction (hyperthyroidism). With the advent of sensitive immunoassays, TSH has become one of the recommended analytes for evaluation of suspected thyroid dysfunction.

Human thyroid stimulating hormone (hTSH) or thyrotropin is a glycoprotein with a molecular weight of approximately 28,000 daltons, synthesized by the basophilic cells (thyrotropes) of the anterior pituitary<sup>2, 8</sup>. hTSH is composed of two non-covalently linked subunits designated alpha ( $\alpha$ ) and beta ( $\beta$ ). The  $\alpha$ -subunit of hTSH is common among several other hormones (follicle stimulating hormone, FSH; luteinizing hormone, LH; and human chorionic gonadotropin, hCG), however, the  $\beta$ -subunit is unique, and confers biological as well as immunological specificity. Both the  $\alpha$ - and  $\beta$ - subunits are required for biological activity<sup>2, 8</sup>.

The FastPack® IPTSH Immunoassay is based on the sandwich immunoassay principle: TSH in the sample (control or patient) is “sandwiched” between the two monoclonal antibodies during the first incubation. The TSH-antibody complex is then bound to the solid-phase paramagnetic particles. After washing the solid-phase to eliminate unbound reagents an alkaline phosphatase substrate is added. The amount of labeled anti-TSH antibody bound to the paramagnetic particles is directly proportional to the concentration of TSH in the sample.

### TEST PRINCIPLE

The FastPack® IPTSH Immunoassay is a chemiluminescence assay based on the “sandwich” principle.

- Primary incubation: Antibody solution (mixture of a biotinylated monoclonal TSH-specific antibody and a monoclonal TSH antibody labeled with alkaline phosphatase) [100  $\mu$ L] reacts with TSH from the patient’s sample, control or calibrator [100  $\mu$ L].
- Secondary incubation: Streptavidin-coated paramagnetic particles are combined with the reaction mixture. During this incubation, the sandwich complex is bound to the solid-phase via the interaction of biotin and streptavidin.
- Removal of unbound materials: The paramagnetic particles are repeatedly washed with wash buffer [0.2 mL/wash] to remove unbound materials.
- Substrate addition and detection: Chemiluminogenic substrate [140  $\mu$ L] is added to the solid-phase bound complex and results in “glow” chemiluminescence, which is measured using the FastPack® IP System and FastPack® System.
- The amount of bound labeled-antibody is directly proportional to the concentration of TSH in the sample.

## REAGENTS – Content and Concentration

Each FastPack® IP carton contains:

- 30 FastPack® IP

Each FastPack® IP TSH Contains:

- Paramagnetic Particles, 150 µL  
Streptavidin-coated paramagnetic particles in buffer containing 0.1% sodium azide as a preservative.
- TSH Antibody Solution, 100 µL  
Antibody solution containing a mixture of a biotinylated mouse monoclonal anti-TSH antibody and a second mouse monoclonal anti-TSH antibody labeled with alkaline phosphatase in a protein matrix containing Proclin® 150 as preservative.
- Wash Buffer, 2.0 mL  
Tris buffer containing surfactants.
- Substrate, 140 µL  
ImmuGlow™: Indoxyl-3-phosphate and lucigenin in buffer containing preservatives.

Materials required but not provided

- FastPack® IP System
- FastPack® TSH Calibrator Kit – Cat. No. 25000024
- FastPack® Control Kit – Cat. No. 25000056(US) or Cat. No. 25000065(International)

## WARNINGS AND PRECAUTIONS

- For *in-vitro* diagnostic use only.
- Do not pipette by mouth.
- Do not eat, drink or smoke in designated work areas.
- Wash hands thoroughly after handling specimen.
- HAMA Interference: some individuals have antibodies to mouse protein (HAMA), which can cause interference in immunoassays that employ antibodies derived from mice<sup>6</sup>.
- FastPack® IP reagents and FastPack® reagents are stable until the expiration date on the label when stored and handled as directed. Do not use FastPack® IP reagents and FastPack® reagents beyond the expiration date.
- Discard used FastPacks into a Biohazard container.
- The components containing sodium azide are classified per applicable European Economic Community (EEC) Directives as: Very toxic and dangerous to the environment (T+ N). The following are appropriate Risk (R) and Safety (S) phrases for sodium azide:

|        |   |
|--------|---|
| R28    | Very toxic if swallowed.  |
| R32    | Contact with acids liberates very toxic gas.  |
| R50/53 | Very toxic to aquatic organisms. May cause long-term adverse effects in the aquatic environment.            |
| S28    | After contact with skin, wash immediately with plenty of soap-suds.   |
| S45    | In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). |
| S60    | This material and its container must be disposed of as hazardous waste.                                     |
| S61    | Avoid release to the environment. Refer to special instructions/safety data sheets.                         |

## STORAGE INSTRUCTIONS

Store at 2 – 8 °C.

## SPECIMEN COLLECTION/PREPARATION

1. Serum and plasma samples can be used for the FastPack® IPTSH Immunoassay.
2. The National Committee for Clinical Laboratory Standards (NCCLS) provides the following recommendations for handling, processing and storing blood.<sup>4, 5</sup>
  - A. Collect all blood samples observing routine precautions for venipuncture.
  - B. For plasma samples:
    - Collect samples in a heparinized (green top) tube.
    - Mix the tube immediately after collection by gently inverting it several times.
    - Plasma should be separated from the cells by centrifugation within 3 hours from time of collection and stored at 2–8 °C. Transfer the plasma from the original tube for storage.
    - If not tested within 24 hours, the sample should be frozen at –20 °C or colder. TSH in plasma, as measured on the FastPack® IP System and Fastpack® System is stable for up to three freeze/thaw cycles.
  - C. For serum samples:
    - Collect samples in serum blood collection tubes.
    - Ensure that complete clot formation has occurred before centrifugation. Some specimens may exhibit increased clotting time, especially those from patients receiving anticoagulant or thrombolytic therapy.
    - Serum should be separated from the clot and centrifuged within 3 hours from time of collection.

- Remove serum from the cells prior to storage at 2–8 °C.
  - If not tested within 24 hours, the sample should be frozen at –20 °C or colder. TSH in serum, as measured on the FastPack® IP System and FastPack® System is stable for up to three freeze/thaw cycles.
- D. Do not freeze samples (–20 °C) for more than two months.
- E. Samples should be free of red blood cells, or other particulate material for optimal results.
- F. Samples showing turbidity and particulate matter should be centrifuged prior to use.
- G. Ensure the samples are free of bubbles.

## ASSAY PROCEDURE

See the FastPack® IP System Procedure Manual for detailed instructions for running the FastPack® IP System and FastPack® System assays.

## INSTRUMENTATION

FastPack® IP System or FastPack® System

## DETAILS OF CALIBRATION

During the FastPack® IP and FastPack® production process, Qualigen generates a master standard curve and places this information in the barcode of each FastPack® IP and FastPack® label, where it can be read by the FastPack® IP System or FastPack® System analyzer during the testing sequence. The analyzer must be calibrated by the end user so that it is properly adjusted for the particular lot of FastPacks that are being used. Separate calibrations must be run for each type of test, i.e. Total PSA, Testosterone, or TSH. The frequency of calibration varies for each test type. For the FastPack® IPTSH Immunoassay, the FastPack® IP System and FastPack® System analyzer must be calibrated once every 14 days or whenever a new lot of TSH FastPacks are to be used.

Whenever the user performs a calibration for a particular lot of FastPacks or uses a new lot of calibrator, 2 FastPacks must be run for calibration (duplicates). When the calibration expires (14 days after initial calibration) 2 FastPacks must be run for calibration. See the FastPack® IP System Procedure Manual for “Running a Calibration”.

Use FastPack® TSH Calibrator Kit – Cat. No. 25000024

## RESULTS

The FastPack® IP System and FastPack® System analyzer uses the information from the barcode to construct a lookup table of x,y values that represent the standard curve and estimates the concentration of unknown samples by linear interpolation.

## QUALITY CONTROL

Quality control materials simulate real specimens and are essential for monitoring the system performance of assays. Good Laboratory Practices (GLP) include the use of control specimens to ensure that all reagents and protocols are performing properly. See FastPack® IP System Procedure Manual for “Control Testing”. At least two levels of quality control materials should be used.

Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.

Controls available: FastPack® Control Kit – Cat. No. 25000056(US) or Cat. No 25000065(International)

## LIMITATION OF PROCEDURE


- Plasma samples to be collected using heparin as the anticoagulant.
- Specimens can be measured within the reportable range of 0.13 µIU/mL and the upper end of the calibration range, 100 µIU/mL.
- Samples >100 µIU/mL should be reported as such or re-run using another method. Samples may be diluted using the zero calibrator (Cal A) until within range of the assay.
- FastPack® IPTSH Immunoassay does not show a high-dose hook effect up to 5,000 µIU/mL.
- Specimen from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits employing mouse monoclonal antibodies.
- Heterophilic antibodies in a sample have the potential to cause interference in immunoassay systems. Infrequently, TSH levels may appear depressed due to heterophilic antibodies present in the patient’s sample or to nonspecific protein binding. If the TSH level is inconsistent with clinical evidence, additional TSH testing is suggested to confirm the result.
- For diagnostic purposes, the FastPack® IPTSH Immunoassay should always be assessed in conjunction with the patient’s medical history, clinical examination and other findings.

- The TSH assay has not been evaluated in Point of Care settings.
- Performance of this assay has not been established with neonatal specimens.

**EXPECTED VALUES**

Samples were obtained from 211 individuals. Samples were obtained from apparently healthy blood donors without any clinically abnormal indications. TSH levels were determined using the FastPack® IPTSH Immunoassay in conjunction with the FastPack® IP System and FastPack® System in order to establish the TSH concentration in the normal population. The normal range (central 95 percentile) for the FastPack® IPTSH Immunoassay is 0.66 – 5.45 µIU/mL. The expected range reflects the donor population of this study group. Each laboratory should determine their own reference range appropriate for their population.

**Frequency Distribution of TSH Values  
Based on Qualigen's FastPack® IP System  
(n = 211)**

**SPECIFIC PERFORMANCE CHARACTERISTICS**

**Precision**

The reproducibility of the TSH assay was measured by testing pooled plasma samples in duplicate for each sample over ten days using three analyzers and two lots of reagents. Precision is expressed as the Coefficient of Variation. The samples used to assess precision were spread out over the linear range of the assay.

| Sample | Number of Samples (N) | Mean (µIU/mL) | Between Run CV% | Between Analyzer CV% | Between Reagent Lot CV% | TOTAL CV% |
|--------|-----------------------|---------------|-----------------|----------------------|-------------------------|-----------|
| 1      | 117                   | 1.54          | 7.4             | 0.5                  | 0.3                     | 9.7       |
| 2      | 117                   | 12.39         | 5.2             | 1.1                  | 4.8                     | 7.5       |
| 3      | 120                   | 46.30         | 5.1             | 5.0                  | 1.1                     | 9.5       |

In addition to the above samples, control materials were tested using three analyzers and two lots of reagents. Precision is expressed as the Coefficient of Variation. The control samples used to assess precision bracket the expected values for apparently health individuals.

| Sample | Number of Samples (N) | Mean (µIU/mL) | Between Run CV% | Between Analyzer CV% | Between Reagent Lot CV% | TOTAL CV% |
|--------|-----------------------|---------------|-----------------|----------------------|-------------------------|-----------|
| 1      | 48                    | 0.54          | 13.8            | 13.3                 | 4.0                     | 20.0      |
| 2      | 48                    | 4.69          | 6.0             | 5.8                  | 4.8                     | 7.0       |

Precision in serum samples was compared to that observed in plasma samples via analysis of matched samples spanning the range of the assay. Precision is expressed as the Coefficient of Variation. The study showed equivalent performance in serum and plasma samples.

| Sample Concentration Range | Number of Samples (N) | Plasma CV | Serum CV |
|----------------------------|-----------------------|-----------|----------|
| 0.0 - 1.0 $\mu$ IU/mL      | 40                    | 6.1%      | 7.0%     |
| 1.0 - 40.0 $\mu$ IU/mL     | 101                   | 4.4%      | 4.4%     |
| 40.0 - 100.0 $\mu$ IU/mL   | 32                    | 5.4%      | 4.0%     |

## RANGE OF LINEARITY

### Serum

Range of linearity in serum was evaluated following the CLSI EP6-A guideline. For TSH in serum as tested by the FastPack<sup>®</sup> IPTSH assay, the method has been demonstrated to be linear from 0.11  $\mu$ IU/mL (limit of detection) to 100  $\mu$ IU/mL within 5  $\mu$ IU/mL in this interval.

### Plasma

Range of linearity in plasma was evaluated following the CLSI EP6-A guideline. For TSH in plasma as tested by the FastPack<sup>®</sup> IPTSH assay, the method has been demonstrated to be linear from 0.08  $\mu$ IU/mL (limit of detection) to 100  $\mu$ IU/mL within 5  $\mu$ IU/mL in this interval.

## METHOD COMPARISONS

### FastPack<sup>®</sup> v. Abbott IMx Ultrasensitive hTSH II


Clinical samples (n= 96 obtained as frozen plasma samples through outside laboratories) were used to compare the plasma values obtained using the FastPack<sup>®</sup> IPTSH method and the plasma values obtained using the Abbott IMx Ultrasensitive hTSH II method. The values were evaluated for agreement using linear regression analysis, with the associated correlation coefficient.

| n  | Range of Observation ( $\mu$ IU/mL) | Intercept ( $\mu$ IU/mL) | Slope | R    |
|----|-------------------------------------|--------------------------|-------|------|
| 96 | 0.03 to 87.99                       | 0.83                     | 1.03  | 0.97 |

### FastPack<sup>®</sup> Serum v. FastPack<sup>®</sup> Plasma

Blood collections were obtained from 74 healthy male and female volunteers between the ages of 18 and 60, and these specimens were processed to heparin plasma and serum samples in parallel. The values were evaluated for agreement using Deming regression analysis, with the associated correlation coefficient.

| n  | Range of Observation ( $\mu$ IU/mL) | Intercept ( $\mu$ IU/mL) | Slope | R    |
|----|-------------------------------------|--------------------------|-------|------|
| 74 | 0.22 – 92.75 (Heparin Plasma)       |                          |       |      |
| 74 | 0.24 – 90.00 (Serum)                | 0.01                     | 0.97  | 1.00 |

### INTERFERING SUBSTANCES

Interfering substances were added to three serum samples and one plasma sample containing known amounts of TSH. The value obtained for the serum or plasma sample with each interfering substance was compared to the value obtained for the serum or plasma without the interfering substance. An acceptance criterion of  $\pm 20\%$  of the value for the serum or plasma sample without added interfering substance was pre-specified for the study. These compounds did not show interference at the levels indicated.

| Serum Samples TSH Concentrations ( $\mu\text{IU/mL}$ ) | Interfering Substance | Concentration, Interfering Substance | % Interference |
|--|-----------------------|--------------------------------------|----------------|
| 1.0, 5.0, 20.0   | Bilirubin             | 40 mg/dL                             | < 10           |
| 1.0, 5.0, 20.0   | Hemoglobin            | 1,000 mg/dL                          | < 20           |
| 1.0, 5.0, 20.0   | Triglycerides         | 1,000 mg/dL                          | < 5            |

| Plasma Sample TSH Concentration ( $\mu\text{IU/mL}$ ) | Interfering Substance | Concentration, Interfering Substance | % Interference |
|---|-----------------------|--------------------------------------|----------------|
| 1.85  | Bilirubin             | 40 mg/dL                             | < 5            |
| 1.63  | Hemoglobin            | 1,000 mg/dL                          | < 5            |
| 1.63  | Triglycerides         | 1,000 mg/dL                          | < 10           |

Interfering substances were added to serum samples containing known amounts of TSH. The value obtained for the serum with each interfering substance was compared to the value obtained for the serum without the interfering substance. These compounds did not show interference at the levels indicated.

| Test Compound | Test Concentration |
|---------------|--------------------|
| d-Biotin      | 250 ng/mL          |

### LIMIT OF BLANK (LOB), LIMIT OF DETECTION (LOD), AND LIMIT OF QUANTITATION (LOQ)

The limit of blank (LOB, the highest measurement likely to be observed for a blank sample), limit of detection (LOD, the lowest amount of analyte in a sample that can be detected with type I and II error rates set to 5%), and limit of quantitation (LOQ, the lowest amount of analyte in a sample that can be reliably detected and at which the total error

meets the pre-specified requirement for accuracy) were determined for serum and plasma according to CLSI EP17-A. In this study, the Limits of Blank for serum and plasma samples were determined from 60 replicate determinations of defibrinated, delipidized, double charcoal filtered serum. Raw RLUs from the assays were converted to apparent  $\mu\text{IU/mL}$  based on the calibration curve for each assay. For each sample type, the LOB was determined as the non-parametric 95<sup>th</sup> percentile of the distribution of values. This value was 0.016  $\mu\text{IU/mL}$  TSH for serum samples and 0.019  $\mu\text{IU/mL}$  TSH for plasma samples.

For serum samples, the LOD was estimated from 60 replicate determinations of three low level serum samples. Based on the Gaussian distributions, LOD was determined by the following equation:

$$\text{LOD} = \text{LOB} + (c_{\beta} \times \text{SD}_S),$$


Where  $c_{\beta} = 1.645/(1-(1/(4 \times f)))$ , where  $f$  is the degrees of freedom, and  $\text{SD}_S$  is the standard deviation of the observations. In this study, the LOD was found to be 0.11  $\mu\text{IU/mL}$  TSH.

For plasma samples, the LOD was estimated from 60 replicate determinations of six low level plasma samples. Based on the non-Gaussian distributions, LOD was determined by identification of the lowest sample with values that displayed <5% overlap with blank values. This concentration was determined to be 0.08  $\mu\text{IU/mL}$  TSH.

For the LOQ determination in serum and plasma, the prospectively defined acceptance criteria for accuracy were 80-120% recovery of target value of low level samples and <20% CV. In this study, a serum sample with an assigned value of 0.127  $\mu\text{IU/mL}$  TSH met the acceptance criteria. For plasma, a sample with an assigned value of 0.129  $\mu\text{IU/mL}$  TSH met the acceptance criteria.


### FUNCTIONAL SENSITIVITY IN SERUM AND HEPARIN PLASMA

Functional sensitivity, defined as the concentration in  $\mu\text{IU/mL}$  TSH at which the assay displays a CV of 20%, was determined both for serum samples and lithium heparin plasma samples. For both methods, the investigation incorporated the preparation of 9 samples in the specified matrix, spanning a range of ~0.05 to 1  $\mu\text{IU/mL}$  TSH, followed by assay of each sample in 30 replicate determinations, (6 replicates/day over 5 days of testing) with one lot of FastPack<sup>®</sup> TSH reagents. Mean, standard deviation (SD), and % CVs were calculated from the replicate determinations. The mean TSH concentration for each sample was then plotted versus % CV and the data fitted via a power function of the form  $Ae^{kx}$ . The concentration at which 20% CV was achieved was then identified from the fitted data. The following chart displays the relationship between mean concentration and % CV. The vertical line in the plot identifies the point at which 20% CV occurs. For the FastPack<sup>®</sup> IPTSH serum assay, this point occurs at 0.08  $\mu\text{IU/mL}$  TSH. Thus, the functional sensitivity of the FastPack<sup>®</sup> IPTSH serum assay is estimated to be 0.08  $\mu\text{IU/mL}$  TSH.

### TSH lithium heparin plasma functional sensitivity

The following chart displays the relationship between mean concentration and % CV. The vertical line in the plot identifies the point at which 20% CV occurs. For the FastPack® IPTSH heparin plasma assay, this point occurs at 0.115  $\mu\text{IU/mL}$  TSH. Thus, the functional sensitivity of the FastPack® IPTSH heparin plasma assay is estimated to be 0.115  $\mu\text{IU/mL}$  TSH.

### SPECIFICITY

#### Serum

Three pools of serum samples were prepared to achieve concentrations of  $\sim 1$ , 5, and 20  $\mu\text{IU/mL}$  TSH. Each sample was separated into four aliquots and three of these were spiked with hCG, LH, and FSH to concentrations of 200,000, 500 and 500 mIU/mL. Samples were tested along with the unspiked reference in five replicates each on two FastPack® instruments using one lot of test packs. Data were analyzed for % difference versus the reference containing no spiked antigens. An acceptance criterion of  $\pm 20\%$  of the unspiked reference value was pre-specified for the study. In no case, do the spiked antigens alter the recovery by more than 20%.

#### Plasma

TSH-stripped plasma was spiked with different concentrations of hCG, FSH, and LH to determine if these materials gave a positive TSH signal in the assay. Cross-reactivity of the antibody used in the FastPack® IPTSH assay was shown to be not detectable for hCG, LH, and FSH tested at 200,000, 500 and 500 mIU/mL, respectively.

## REFERENCES

- <sup>1</sup> Braverman, LE, Utigen, RD, Eds: Werner and Ingbar's The Thyroid. A Fundamental and Clinical Text. 7th ed., Philadelphia, Lippincott-Raven, 1996.
- <sup>2</sup> Burtis, CA, Ashwood, ER, Eds: Tietz Textbook of Clinical Chemistry, 3rd ed., New York, WB Saunders Company, p. 1492, 1999.
- <sup>3</sup> Larsen, PR: Thyroid-pituitary interaction. N. Engl. J. Med., 306:23-32, 1982.
- <sup>4</sup> Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5<sup>th</sup> Edition: H3-A5: 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS)
- <sup>5</sup> Approved guideline – procedures for the handling and processing of blood specimens, H18-A2; 19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- <sup>6</sup> Schroff, RJ, Foon, KA, *et.al*: Human anti-mouse immunoglobulin responses in patients receiving monoclonal antibody therapy. Cancer Res, 45:879 – 885, 1985
- <sup>7</sup> Jackson, IM: Thyrotropin-releasing hormone. N. Engl. J. Med., 306: 145-155, 1982.
- <sup>8</sup> Pierce, JG: Subunits of Pituitary Thyrotropin. Their relationship to other glycoprotein hormones. Endocrinology, 89: 1331-1344, 1971.
- <sup>9</sup> Nicoloff JT, Spencer CA: The use and misuse of the sensitivity thyrotropin assays. J. Clin. Endocrinology, 71, 553 – 558, 1990.
- <sup>10</sup> Spencer CA: Thyroid profiling for the 1990's: free T4 estimate or sensitive TSH measurement. J. Clin Immunoassay, 12, 82 – 89, 1989.

Qualigen, Inc.  
Carlsbad, CA 92011 USA  
Technical Support  
(760) 918-9165  
(877) 709-2169

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© 2000 Qualigen, Inc. All rights reserved. Qualigen and FastPack are trademarks or registered trademarks of Qualigen, Inc. All other trademarks are the property of their respective owners.

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## REFERENCIAS

- <sup>1</sup> Braverman, LE, Utigen, RD, Eds: Werner and Ingbar's The Thyroid. A Fundamental and Clinical Text. 7th ed., Philadelphia, Lippincott-Raven, 1996.
- <sup>2</sup> Burtis, CA, Ashwood, ER, Eds: Tietz Textbook of Clinical Chemistry, 3rd ed., New York, WB Saunders Company, p. 1492, 1999.
- <sup>3</sup> Larsen, PR: Thyroid-pituitary interaction. N. Engl. J. Med., 306:23-32, 1982.
- <sup>4</sup> Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5<sup>th</sup> Edition: H3-A5: 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS)
- <sup>5</sup> Approved guideline – procedures for the handling and processing of blood specimens, H18-A2; 19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- <sup>6</sup> Schrott, RJ, Foon, KA, *et.al.*: Human anti-mouse immunoglobulin responses in patients receiving monoclonal antibody therapy. Cancer Res, 45:879 – 885, 1985
- <sup>7</sup> Jackson, IM: Thyrotropin-releasing hormone. N. Engl. J. Med., 306: 145-155, 1982.
- <sup>8</sup> Pierce, JG: Subunits of Pituitary Thyrotropin. Their relationship to other glycoprotein hormones. Endocrinology, 89: 1331-1344, 1971.
- <sup>9</sup> Nicoloff JT, Spencer CA: The use and misuse of the sensitivity thyrotropin assays. J. Clin. Endocrinology, 71, 553 – 558, 1990.
- <sup>10</sup> Spencer CA: Thyroid profiling for the 1990's: free T4 estimate or sensitive TSH measurement. J. Clin Immunoassay, 12, 82 – 89, 1989.

Qualigen, Inc.  
Carlsbad, CA 92011 EE.UU.  
Soporte Técnico  
+1(760) 918-9165  
+1(877) 709-2169

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