CAUTION: United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.

INTENDED USE

The FastPack® IP hCG Immunoassay is a paramagnetic particle immunoassay for the *in-vitro* quantitative determination of hCG in human serum or EDTA or lithium-heparin plasma. The FastPack® IP hCG Immunoassay is designed for use with the FastPack® IP System and is indicated for the early detection of pregnancy.

SUMMARY

Human chorionic gonadotropin (hCG) is an oligosaccharide glycoprotein hormone produced by the placenta. It is a heterodimeric protein with an alpha subunit chain of 92 amino acids in length with a molecular weight of 15,000-20,000 Daltons and a beta subunit chain of 152 amino acids in length with a molecular weight of 25,000-30,000 Daltons. The amino acid sequence of the alpha subunit is identical to three hormones produced by the pituitary gland: luteinizing hormone (LH), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH). However, the amino acid sequence of the beta subunit of hCG is unique, as are the beta subunits of LH, FSH, and TSH, and thus confers immunological specificity to assays developed to measure the four hormones.

Shortly after implantation of a fertilized ovum into the uterine wall, the trophoblast begins to produce hCG and continues to do so until the placenta is capable of production and secretion of hCG¹. During these initial stages of pregnancy, hCG promotes the maintenance of the corpus luteum by stimulation of progesterone secretion. During a normal pregnancy, hCG measured in serum varies in concentration but may be as high as 50 mIU/mL in the week following conception, and doubles every 1.5-3 days for the subsequent six weeks². Serum hCG concentrations continue to rise until the end of the first trimester, after which concentrations diminish for the remainder of the gestation period. Following delivery, hCG concentrations in serum rapidly return to pre-pregnancy levels.

Due to the very rapid rise in hCG production in the first trimester, hCG has been historically measured as an indicator of pregnancy in humans. The elevation of hCG levels is quite specific for pregnancy as healthy, non-pregnant individuals have low to undetectable concentrations of hCG in serum. In a confirmed pregnancy, unusually low levels of hCG may also result from other conditions such as ectopic pregnancy or impending spontaneous abortions³.

hCG was originally measured by correlation to gonadal tissue response in various animals. Later, tests to measure hCG in human urine collections were developed, typically employing latex agglutination methodologies. Radioimmunoassay (RIA) techniques based on competitive inhibition methods utilizing polyclonal antibodies were developed in the early 1970's and were shown to be sensitive for measurements of hCG in serum and urine⁴. Subsequently, two-site immunoradiometric (IRMA) assays employing dual monoclonal antibodies specific for hCG were developed and provided a platform for measurement of hCG with increased sensitivity, specificity, and precision⁵.

TEST PRINCIPLE

The FastPack® IP hCG Immunoassay is a chemiluminescence assay based on the "sandwich" principle.

- Primary incubation: A monoclonal anti-hCG antibody labeled with alkaline phosphatase [100 µL] reacts with hCG from the patient's sample, control or calibrator [25 µL].
- Secondary incubation: A second monoclonal anti-hCG antibody covalently coupled to biotin and pre-bound to streptavidin-coated paramagnetic particles [150 µL] is combined with the reaction mixture.
- Removal of unbound materials: The paramagnetic particles are repeatedly washed with wash buffer [0.2 mL/wash] to remove unbound materials.
- Substrate addition and detection: Chemiluminogenic substrate [140 µL] is added to the solid-phase bound complex and results in "glow" chemiluminescence, which is measured using the FastPack® IP System.
- The amount of bound labeled-antibody is directly proportional to the concentration of hCG in the sample.

REAGENTS – Content and Concentration

Each FastPack® IP carton contains:

- 30 FastPack® IPs

Each FastPack® IP Contains:

- Paramagnetic Particles, 150 µL
Biotin-mouse monoclonal anti-βhCG antibody bound to streptavidin-coated paramagnetic particles in buffer containing <0.1% sodium azide and 0.1% ProClin® 300 as a preservative.
- hCG Antibody Solution 1, 100 µL
Antibody solution containing a mouse monoclonal anti-βhCG antibody labeled with alkaline phosphatase in a protein matrix containing <0.1% sodium azide and 0.1% ProClin® 300 as a preservative.
- Wash Buffer, 2.0 mL
Tris buffer containing surfactants.
- Substrate, 140 µL
ImmuGlow™ Plus: Indoxyl-3-phosphate and lucigenin in buffer containing preservatives.

Materials required but not provided

- FastPack® IP System
- FastPack® hCG Calibrator Kit – Cat. No. 25000035
- FastPack® hCG Sample Diluent – Cat. No. 25000036
- Quality control materials (see “Quality Control” section)

WARNINGS AND PRECAUTIONS

- For *in-vitro* diagnostic use only.
- Do not pipette by mouth.
- Do not eat, drink or smoke in designated work areas.
- Wash hands thoroughly after handling specimen.
- HAMA Interference: some individuals have antibodies to mouse protein (HAMA), which can cause interference in immunoassays that employ antibodies derived from mice⁶.
- FastPack® IP reagents are stable until the expiration date on the label when stored and handled as directed. Do not use FastPack® IP reagents beyond the expiration date.
- Discard used FastPacks into a Biohazard container.
- ProClin® 300 is an irritant. The following are appropriate Risk (R) and Safety (S) phrases for ProClin® 300:
 - R43 May cause sensitization by skin contact
 - S28-37 After contact with skin, wash immediately with plenty of soap and water. Wear suitable gloves.
- The components containing sodium azide are classified per applicable European Economic Community (EEC) Directives as: Very toxic and dangerous to the environment (T+ N). The following are appropriate Risk (R) and Safety (S) phrases for sodium azide:
 - R28 Very toxic if swallowed.
 - R32 Contact with acids liberates very toxic gas.
 - R50/53 Very toxic to aquatic organisms. May cause long-term adverse effects in the aquatic environment.
 - S28 After contact with skin, wash immediately with plenty of soap-suds.
 - S45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
 - S60 This material and its container must be disposed of as hazardous waste.
 - S61 Avoid release to the environment. Refer to special instructions/safety data sheets.

STORAGE INSTRUCTIONS

Store at 2 – 8 °C.

SPECIMEN COLLECTION/PREPARATION

1. Serum, EDTA or lithium-heparin plasma samples can be used for the FastPack® IP hCG Immunoassay.
2. The National Committee for Clinical Laboratory Standards (NCCLS) provides the following recommendations for handling, processing and storing blood.^{7,8}
 - A. Collect all blood samples observing routine precautions for venipuncture.
 - B. For plasma samples:
 - Collect samples in an EDTA (lavender top) or heparinized (green top) tube.
 - Mix the tube immediately after collection by gently inverting it several times.
 - Plasma should be separated from the cells by centrifugation within 3 hours from time of collection and stored at 2–8° C. Transfer the plasma from the original tube for storage.
 - If not tested within 24 hours, the sample should be frozen at –20° C or colder.

- C Samples should be free of red blood cells, or other particulate material for optimal results.
- D Samples showing turbidity and particulate matter should be centrifuged prior to use.
- E Ensure the samples are free of bubbles.

ASSAY PROCEDURE

See the FastPack® IP System Procedure Manual for detailed instructions for running FastPack® IP assays.

INSTRUMENTATION

FastPack® IP System

DETAILS OF CALIBRATION

During the FastPack® IP and FastPack® production process, Qualigen generates a master standard curve and places this information in the barcode of each FastPack® IP FastPack® label, where it can be read by the FastPack® IP System analyzer during the testing sequence. The FastPack® IP System analyzer must be calibrated by the end user so that it is properly adjusted for the particular lot of FastPacks that are being used. Separate calibrations must be run for each type of test, i.e. Total PSA, Testosterone, or hCG. The frequency of calibration varies for each test type. For the FastPack® IP hCG Immunoassays, the FastPack® IP System analyzer must be calibrated once every 30 days or whenever a new lot of hCG FastPacks is to be used.

Whenever the user performs a calibration for a particular lot of FastPacks or uses a new lot of calibrator, 2 FastPacks must be run for calibration (duplicates). When the calibration expires (30 days after initial calibration) 2 FastPacks must be run for calibration. See FastPack® IP System Procedure Manual for "Running a Calibration".

Use FastPack® hCG Calibrator Kit – Cat. No. 25000035

RESULTS

The FastPack® IP System analyzer uses the information from the barcode to construct a lookup table of x,y values that represent the standard curve and estimates the concentration of unknown samples by linear interpolation.

QUALITY CONTROL


Quality control materials are essential for monitoring the system performance of assays. Good Laboratory Practices (GLP) include the use of control specimens to ensure that all reagents and protocols are performing properly. See FastPack® IP System Procedure Manual for "Control Testing". At least two levels of quality control materials should be used with hCG concentrations at or near medical decision levels. Please contact Qualigen Customer Service for recommendations of quality control materials. Users should also follow the appropriate federal, state and local guidelines concerning the running of external quality controls.

LIMITATIONS OF PROCEDURE

- Plasma samples to be collected using EDTA or lithium-heparin as the anticoagulants.
- Specimens can be detected within the reportable range of the limit of detection (1.8 mIU/mL hCG) and the upper end of the calibration range, 1000 mIU/mL and reliably measured at concentrations at and above the LOQ (3.1 mIU/mL). Values below the LOD should be reported as "<1.8 mIU/mL hCG" and values above the range should be reported as ">1000 mIU/mL hCG" if not diluted and re-tested.
- Samples >1000 mIU/mL should be reported as such or re-run using another method. Samples may be diluted using the hCG sample diluent until within range of the assay.
- FastPack® IP hCG Immunoassay does not show a high-dose hook effect up to 500,000 mIU/mL.
- Specimen from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits employing mouse monoclonal antibodies.
- Heterophilic antibodies in a sample have the potential to cause interference in immunoassay systems. Infrequently, hCG levels may appear depressed due to heterophilic antibodies present in the patient's sample or to nonspecific protein binding. If the hCG level is inconsistent with clinical evidence, additional hCG testing is suggested to confirm the result.
- For diagnostic purposes, the FastPack® IP hCG Immunoassay should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.
- The hCG assay has not been evaluated in Point of Care settings.
- The concentration of Human Chorionic Gonadotropin (hCG) in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the hCG assay method used. Values obtained with different assay methods should not be used interchangeably.

EXPECTED VALUES/REFERENCE INTERVALS

Samples were obtained from 241 individuals, 115 males and 126 non-pregnant females. Samples were obtained from apparently healthy blood donors without any clinically abnormal indications. hCG levels were determined using the FastPack® IP hCG Immunoassay in conjunction with the FastPack® IP System in order to establish the hCG concentrations in the reference populations. The reference interval (5th to 95th percentiles) for the FastPack® IP hCG Immunoassay for males is <1.8 mIU/mL and <1.8 – 3.0 mIU/mL for non-pregnant females. These reference intervals reflect the donor population of this study group. Each laboratory should determine their own reference intervals appropriate for their population. The distribution of values for non-pregnant females is presented in the histogram below.

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

Three control serum pools spiked with varying amounts of hCG were assayed using two lots of FastPack® IP hCG reagents and three FastPack® IP instruments. Each control was tested in duplicate determinations per day over a period of 10 days for a total of 20 replicate determinations per control. The data was analyzed using the analysis of variance (ANOVA) technique and within-assay and total imprecision were determined. Mean, standard deviation (SD), and percent coefficient of variation (%CV) were calculated for each control solution.

Within-Assay Imprecision

	Control 1	Control 2	Control 3	Control 4
Mean mIU/mL hCG	3.4	19.0	347.3	722.3
SD	0.37	1.8	21.5	73.7
%CV	10.8	9.4	6.2	10.2

Total Imprecision

	Control 1	Control 2	Control 3	Control 4
Mean mIU/mL hCG	3.4	19.0	347.3	722.3
SD	0.50	2.2	27.1	98.3
%CV	14.5	11.4	7.8	13.6

Range of linearity

For hCG as tested by the FastPack® IP hCG assay, based on the CLSI EP-6A evaluation the method has been demonstrated to be linear from 1.8 mIU/mL to 940.6 mIU/mL within 5 mIU/mL in the interval between 1.8 and 50 mIU/mL and 25 mIU/mL in the interval between 50 and 1000 mIU/mL hCG.

Linearity Upon Dilution

Two serum pools with high endogenous concentrations of hCG were diluted 1:2, 1:4, 1:8, and 1:16 with a protein-based Tris buffer diluent and quantified in triplicate determinations with one lot of FastPack® IP hCG reagents on a single FastPack® IP instrument. Measured values of the dilutions were compared to their expected values (determined in the same assay) based on the measured value of the Neat sample and the respective dilution factors. (% Recovery = [Observed value/Expected Value] X 100).


Serum Pool 1	Expected (mIU/mL)	Measured (mIU/mL)	% Recovery
Neat	N/A	699.2	N/A
1:2	349.6	366.6	104.9
1:4	174.8	174.7	99.9
1:8	87.4	90.1	103.1
1:16	43.7	46.5	106.4

Serum Pool 2	Expected (mIU/mL)	Measured (mIU/mL)	% Recovery
Neat	N/A	443.1	N/A
1:2	221.5	205.5	92.7
1:4	110.8	105.3	95.0
1:8	55.4	59.2	107.0
1:16	27.7	28.7	103.7

Method Comparison

Clinical serum samples (n=106) were used to compare the values obtained using the FastPack® IP hCG method and the values obtained using the Beckman-Coulter Access® hCG method. The values were evaluated for agreement using Deming regression analysis, with the associated correlation coefficient.

n	Range of Observation (mIU/mL)	Intercept (mIU/mL)	Slope	R ²
106	0.0 – 864.6	2.44	1.11	0.99

INTERFERING SUBSTANCES

The effect of icteric, lipemic, and hemolyzed serum samples and serum samples with high protein content on quantification of hCG was investigated by spiking of one serum pool with endogenous hCG concentration of ~300 mIU/mL with known concentrations of bilirubin, triglycerides, hemoglobin, and human serum albumin and another serum pool with endogenous hCG concentration of ~3 mIU/mL with the same concentrations of the interfering substances. The values obtained for the serum pools with each interfering substance were compared to the values obtained for the serum pools without the interfering substance and the percentage bias or absolute bias in mIU/mL determined. These compounds did not show interference at the levels indicated.

	Interfering Substance			
	Bilirubin (30 mg/dL)	Lipemia (1800 mg/dL)	Hemoglobin (500 mg/dL)	Human Albumin (3 g/dL)
Non-spiked aliquot	303.4 mIU/mL	274.6 mIU/mL	260.6 mIU/mL	255.0 mIU/mL
Spiked aliquot	318.0 mIU/mL	288.9 mIU/mL	265.9 mIU/mL	260.6 mIU/mL
% Bias	+4.8%	+5.2%	+5.3%	+2.2%

	Interfering Substance			
	Bilirubin (30 mg/dL)	Lipemia (1800 mg/dL)	Hemoglobin (500 mg/dL)	Human Albumin (3 g/dL)
Non-spiked aliquot	3.8 mIU/mL	3.3 mIU/mL	2.3 mIU/mL	4.3 mIU/mL
Spiked aliquot	4.2 mIU/mL	3.7 mIU/mL	1.9 mIU/mL	4.4 mIU/mL
Absolute bias	0.4	0.4	-0.4	0.1

Interfering substances were added to serum samples containing known amounts of hCG. The value obtained for the serum with each interfering substance was compared to the value obtained for the serum without the interfering substance. These compounds did not show interference at the levels indicated.

Test Compound	Test Concentration
Biotin	2000 ng/mL

Cross-reactivity

No significant cross reactivity (<0.5 mIU/mL) was obtained when FastPack® hCG sample diluent was spiked with the following cross-reactants: LH, at 1000 mIU/mL, FSH at 1000 mIU/mL and TSH at 1000 µIU/mL. The molar percent specificity when 25 mIU/mL of WHO 75/551 free βhCG subunit is added to FastPack® hCG sample diluent is approximately 110%.

Limit of Blank (LOB), Limit of Detection (LOD), and Limit of Quantitation (LOQ)

The limit of blank (LOB, the highest measurement likely to be observed for a blank sample), limit of detection (LOD, the lowest amount of analyte in a sample that can be detected with type I and II error rates set to 5%), and limit of quantitation (LOQ, the lowest amount of analyte in a sample that can be reliably detected and at which the total error meets the pre-specified requirement for accuracy) were determined according to CLSI EP17-A for the FastPack® IP hCG assay. In this study, the limit of blank was determined from 10 replicate determinations of the FastPack® hCG sample diluent on each of three different FastPack® IP instruments using two different lots of reagent sets. Raw RLU values from the assays were converted to apparent mIU/mL based on the calibration curve for each assay. The LOB was determined as the maximum observed value. This value was 0.8 mIU/mL hCG.

The LOD was estimated from 80 replicate determinations of a low control sample. Per the CLSI EP17-A guideline, LOD was determined by the following equation:

$$\text{LOD} = \text{LOB} + (c_{\beta} * \text{SDs}),$$

Where $c_{\beta} = 1.645/(1-(1/(4 * f)))$, where f is the degrees of freedom, and SDs is the standard deviation of the observations. In this study, the LOD was found to be 1.8 mIU/mL hCG.

For the LOQ analyses, the prospectively defined goals for accuracy were 80-120% recovery and <20% CV. In this study, a control sample recovering at 3.1 mIU/mL hCG was shown to meet these criteria for recovery and CV (93.8% and 18.7%, respectively), thus the LOQ was found to be 3.1 mIU/mL hCG.

REFERENCES

- ¹ Kardana A, et al. The heterogeneity of hCG. *Endocrinology* 1991; 129: 1541-1567.
- ² Vaitukaitis JL. Recent progress in hormone research 1976; 32: 289.
- ³ Sokolove PJ, Faix JD. Agreement of intact and beta chain-specific HCG assays in abnormal pregnancy. *Journal of Clinical Immunoassay*, Fall, 1991; 14, No. 3: 196-199.
- ⁴ Vaitukaitis JL, Braunstein GD, Ross GT. A radioimmunoassay which specifically measures human chorionic gonadotropin in the presence of human luteinizing hormone. *American Journal of Obstetrics and Gynecology* 1972; 113: 751-758.
- ⁵ Norman RJ, Buck RH, Joubert SM. Comparison of human chorionic gonadotrophin concentrations in the sera of patients with normal and abnormal pregnancy measured by radioimmunoassay and immunoradiometric assay. *South African Medical Journal* April, 1989; 75: 318-319.
- ⁶ Schroff, RJ, Foon, KA, *et.al*: Human anti-mouse immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res*, 45:879 – 885, 1985
- ⁷ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition: H3-A5: 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS)
- ⁸ Approved guideline – procedures for the handling and processing of blood specimens, H18-A2; 19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).

Qualigen, Inc.
Carlsbad, CA 92011 USA
Technical Support
(760) 918-9165
(877) 709-2169

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© 2000 Qualigen, Inc. All rights reserved. Qualigen and FastPack are trademarks or registered trademarks of Qualigen, Inc.
All other trademarks are the property of their respective owners.

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REFERENZEN

- ¹ Kardana A, et al. The heterogeneity of hCG. *Endocrinology* 1991; 129: 1541-1567.
- ² Vaitukaitis JL. Recent progress in hormone research 1976; 32: 289.
- ³ Sokolove PJ, Faix JD. Agreement of intact and beta chain-specific HCG assays in abnormal pregnancy. *Journal of Clinical Immunoassay*, Fall, 1991; 14, No. 3: 196-199.
- ⁴ Vaitukaitis JL, Braunstein GD, Ross GT. A radioimmunoassay which specifically measures human chorionic gonadotropin in the presence of human luteinizing hormone. *American Journal of Obstetrics and Gynecology* 1972; 113: 751-758.
- ⁵ Norman RJ, Buck RH, Joubert SM. Comparison of human chorionic gonadotrophin concentrations in the sera of patients with normal and abnormal pregnancy measured by radioimmunoassay and immunoradiometric assay. *South African Medical Journal* April, 1989; 75: 318-319.
- ⁶ Schroff, RJ, Foon, KA, *et.al.*: Human anti-mouse immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res*, 45:879 – 885, 1985
- ⁷ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition: H3-A5: 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS)
- ⁸ Approved guideline – procedures for the handling and processing of blood specimens, H18-A2; 19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).

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REFERENCIAS

- ¹ Kardana A, et al. The heterogeneity of hCG. *Endocrinology* 1991; 129: 1541-1567.
- ² Vaitukaitis JL. Recent progress in hormone research 1976; 32: 289.
- ³ Sokolove PJ, Faix JD. Agreement of intact and beta chain-specific HCG assays in abnormal pregnancy. *Journal of Clinical Immunoassay*, Fall, 1991; 14, No. 3: 196-199.
- ⁴ Vaitukaitis JL, Braunstein GD, Ross GT. A radioimmunoassay which specifically measures human chorionic gonadotropin in the presence of human luteinizing hormone. *American Journal of Obstetrics and Gynecology* 1972; 113: 751-758.
- ⁵ Norman RJ, Buck RH, Joubert SM. Comparison of human chorionic gonadotrophin concentrations in the sera of patients with normal and abnormal pregnancy measured by radioimmunoassay and immunoradiometric assay. *South African Medical Journal* April, 1989; 75: 318-319.
- ⁶ Schroff, RJ, Foon, KA, *et.al.* Human anti-mouse immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res*, 45:879 – 885, 1985
- ⁷ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition: H3-A5: 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS)
- ⁸ Approved guideline – procedures for the handling and processing of blood specimens, H18-A2; 19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).

Qualigen, Inc.
Carlsbad, CA 92011 EE.UU.
Soporte Técnico:
+1 (760) 918-9165
+1 (877) 709-2169

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