For the Quantitative Measurement of Prostate-Specific Antigen (PSA) in Human Serum and Plasma

The concentration of PSA in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the PSA assay method used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining PSA levels serially is changed, additional sequential testing should be carried out to confirm baseline values.

PSA concentrations are dependent on the standard used to calibrate the assay. PSA concentrations based on calibration to the WHO 96/670 Reference Preparation will differ significantly from PSA concentrations based on calibration to the original Hybritech assay. The concentrations are not inter-changeable. If the calibration is changed, accepted laboratory practice is to establish a new baseline for patient monitoring.¹ For more information on WHO and Hybritech configuration reference your instrument procedure manual.

INTENDED USE

The FastPack® IP Total PSA Immunoassay is a paramagnetic particle immunoassay for the *in vitro* quantitative determination of prostate-specific antigen (PSA) in human serum and plasma as an aid in the management of patients with prostate cancer. The FastPack® IP Total PSA Immunoassay is designed for use with the FastPack® IP System.

SUMMARY

Prostate-specific antigen (PSA) is a glycoprotein (molecular weight 30,000-34,000 Daltons), which exhibits a high degree of homology with serine proteases of the kallikrein family. It has the function of serine protease with chymotrypsin-like activity.² The proteolytic activity of PSA in blood is inhibited by the irreversible formation of complexes with protease inhibitors such as α -1- antichymotrypsin (ACT), α -2-macroglobulin and other acute phase proteins.³ PSA occurs in three forms in blood. The two immunodetectable forms are PSA complexed with the serine protease inhibitor α -1-antichymotrypsin and free or uncomplexed PSA.^{4, 5, 6} The third form, which is undetectable in conventional immunoassay, is PSA complexed with α -2 macroglobulin.

This is due to engulfment and masking of PSA epitopes by the α -2-macroglobulin molecule. Elevated levels of PSA in serum are generally indicative of a pathological condition of the prostate, for example, carcinoma, benign prostatic hyperplasia (BPH) or prostatitis.^{7,8} Prostate Specific Antigen is also present in paraurethral and anal glands, as well as breast tissue. An inflammation, trauma or stimulation of the prostate (e.g. following biopsy and colonoscopy, etc.) can result in PSA elevations of varying duration and magnitude.

Determination of PSA is an integral part in monitoring the efficiency of therapy in patients with prostate cancer, for example, radical prostatectomy, radiotherapy or patients receiving hormonal therapy.⁹

TEST PRINCIPAL

The FastPack® IP Total PSA Immunoassay is a chemiluminescence assays based on the “sandwich” principle.

- Primary incubation: Specimen, control or calibrator [25 μ L] and antibody solution (mixture of a biotinylated monoclonal PSA-specific antibody and a monoclonal PSA-specific antibody labeled with alkaline phosphatase) [100 μ L] react to form a sandwich complex.
- Secondary incubation: Streptavidin-coated paramagnetic particle solution is added to the reaction mixture. During this incubation, the sandwich complex is bound to the solid-phase via the interaction of biotin and streptavidin.
- Removal of unbound materials: The paramagnetic particles are washed with wash buffer [0.2 mL/wash] to remove unbound materials.
- Substrate addition and detection: Chemiluminogenic substrate [140 μ L] is added to the solid-phase bound complex and results in “glow” chemiluminescence, which is measured using the FastPack® IP System analyzer.
- The amount of bound labeled-antibody is directly proportional to the concentration of PSA in the sample.

REAGENTS - Content and Concentration

FastPack® IP Total PSA Kit - Cat. No. 25000079

Each FastPack® IP kit contains:

- 30 FastPacks
- 1 Control 1 vial x 5mL
- 1 Control 2 vial x 5mL
- 2 Control Range Cards
 - For TPSA the controls and ranges are to be used only with the TPSA FastPacks in which they are kitted
 - For FPSA the controls and ranges may be use with any lot of FPSA

Each FastPack® IP Contains:

- Paramagnetic Particles, 150 µL
Streptavidin-coated paramagnetic particles in buffer containing 0.1% sodium azide as a preservative.
- PSA Antibody Solution, 100 µL
Antibody solution containing mouse monoclonal antibody coupled to biotin and mouse monoclonal antibody labeled with alkaline phosphatase in a protein matrix containing 0.1% sodium azide as preservative.
- Wash Buffer, 2.0 mL
Tris buffer containing surfactants.
- Substrate, 145 µL
ImmuGlow™: Indoxyl-3-phosphate and lucigenin in buffer containing preservatives.

Control Information:

- 5mL / vial. Liquid.
- Components of human origin prepared in a Tris Buffer with protein stabilizers and 0.1% Sodium Azide.
- The total PSA target values and ranges are printed on the Control Range Card.
- Mix contents by gently inverting before use. Avoid bubble formation.

Materials required but not provided

- FastPack® IP System
- FastPack® Total PSA Calibrator Kit - Cat. No. 25000077

WARNING AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- **Human source material. Treat as potentially infectious.**
- Do not eat, drink or smoke in designated work areas.
- Do not mix controls from different lots.
- After opening, controls remain stable until the expiration date on the label when stored and handled as directed. Do not use controls beyond the expiration date.
- Avoid microbial contamination of reagent when removing aliquots from the bottles.
- Wash hands thoroughly after handling specimen.
- HAMA Interference: some individuals have antibodies to mouse protein (HAMA), which can cause interference in immunoassays that employ antibodies derived from mice. In particular, it has been reported that serum or plasma samples from patients who have undergone therapy or diagnostic procedures that include infusion of mouse monoclonal antibody may produce erroneous results in such assays.
- FastPack® IP reagents are stable until the expiration date on the label when stored and handled as directed. Do not use FastPack® IP reagents beyond expiration date.
- Discard used FastPacks into a biohazard container.
- Discard unused or expired control material, in stoppered vial, into a Biohazard container. The components containing sodium azide are classified per applicable European Economic Community (EEC) Directives as: Harmful (Xn). The following are appropriate Risk (R) and Safety (S) phrases:

R28	Very Toxic if swallowed.
R32	Contact with acids liberates very toxic gas.
R50/53	Very Toxic to aquatic organisms may cause long-term adverse effects in the aquatic environment.
S2	Keep out of reach of children.
S13	Keep away from food, drink and animal feeding stuffs.
S36	Wear suitable protective clothing.
S46	If swallowed, seek medical advice immediately and show this container or label.

STORAGE INSTRUCTIONS

Store at 2 - 8 °C.

SPECIMEN COLLECTION/PREPARATION

1. Serum or plasma samples can be used for the FastPack® IP Total PSA Immunoassay.
2. The National Committee for Clinical Laboratory Standards (NCCLS) provides recommendations for handling, processing, and storing blood.^{10,11}
3. Collect all blood samples observing routine precautions for venipuncture.
4. It is not required that patients fast prior to blood collection.
5. For serum samples:
 - Ensure that complete clot formation has occurred before centrifugation. This takes approximately 30 minutes. Some samples may exhibit increased clotting time, especially those from patients receiving anticoagulant or thrombolytic therapy.
 - Serum should be centrifuged and separated from the clot within 3 hours from time of collection.
 - Remove serum from the cells prior to storage at 2-8 °C.
 - If not tested within 24 hours, the sample should be frozen at –20 °C or colder.
6. For plasma samples:
 - Collect samples in Lithium Heparin tubes. Citrate tubes should not be used.
 - Plasma should be centrifuged and separated within 3 hours from time of collection.
 - Remove plasma from the cells prior to storage at 2-8 °C.
 - If not tested within 24 hours, the sample should be frozen at –20 °C or colder.
7. Do not freeze samples (-20° C) for more than two months.
8. Frozen samples should be thawed completely and mixed by gentle inversion prior to use.
9. Samples should be free of Fibrin, Red Blood Cells, or other particulate material for optimal results. Samples showing turbidity and/or particulate matter should be centrifuged prior to use.
10. Ensure the samples are free of bubbles.
11. Samples collected up to 2 hours after DRE show no significant increases in PSA.¹⁹
12. Human samples should be handled in accordance with the OSHA standard on Bloodborne Pathogens.²⁰

ASSAY PROCEDURE

See the FastPack® IP System Procedure Manual for information on operating the FastPack® IP System.

INSTRUMENTATION

FastPack® IP System

DETAILS OF CALIBRATION

During the FastPack® IP production process, Qualigen generates a master standard curve and places this information in the barcode of each FastPack® IP label, where it can be read by the FastPack® IP System analyzer during the testing sequence. The FastPack® IP System analyzer must be calibrated by the user to ensure that it is properly adjusted for the particular lot of FastPacks that is being used. Separate calibrations must be run for each type of test, i.e. Free PSA, Total PSA or Testosterone. The frequency of calibration varies for each test type. For the FastPack® IP Total PSA Immunoassay, the FastPack® IP System analyzer must be calibrated once every 30 days or whenever a new lot of Total PSA FastPacks are to be used.

Whenever the user performs an initial calibration for a particular lot of FastPacks or uses a new lot of calibrator, 2 FastPacks must be run for calibration (duplicates). Whenever recalibration is performed with the same lot of FastPacks and calibrator, 2 FastPacks must be run for calibration. See FastPack® IP System Procedure Manual for “Calibrating the FastPack® IP System”.

Use FastPack® Total PSA Calibrator Kit – Cat. No. 25000077

RESULTS

The FastPack® IP System analyzer uses the information from the barcode to construct a lookup table of x,y values that represent the standard curve and estimates the concentration of unknown samples by linear interpolation.

QUALITY CONTROL

Quality control materials simulate real specimens and are essential for monitoring the system performance of assays. Good Laboratory Practices (GLP) include the use of control specimens to ensure that all reagents and protocols are performing properly. See FastPack® IP System Procedure Manual for “Control Testing”.

The included PSA controls are assayed quality control materials used for the verification of the accuracy and precision of the FastPack® IP System when used for the quantitative determination of PSA in human serum and plasma. The use of control material is indicated as an objective assessment of the precision of methods and techniques in use. Two levels of control are provided to allow performance monitoring within the clinical range.

LIMITATION OF PROCEDURE

- The FastPack® IP Total PSA is not intended to be used to diagnose prostate cancer.
- Samples can be accurately measured within the reportable range of the analytical sensitivity and the highest calibrator, 50 ng/mL (40 ng/mL WHO).
- Samples >50 ng/mL (>40 ng/mL WHO), or where there is the clinical possibility of a “high-dose hook effect”, should be run using another method. Dilution of out of range results is not recommended.
- In a two-site immunoassay, samples with extremely high concentrations may paradoxically read within the calibration range of the assay. This possibility should be considered in patients with disseminated cancer of the prostate manifesting clinically inappropriately low total PSA concentrations. The FastPack® IP Total PSA Immunoassay does not show a high dose hook effect up to 7,500 ng/mL.
- Samples from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such samples may show either falsely elevated or depressed values when tested with assay kits employing mouse monoclonal antibodies.^{12,13}
- It is known that in rare cases PSA isoforms do exist which may be measured differently by different PSA tests. Findings of this kind have occasionally been reported for PSA tests from various manufacturers.^{14, 15, 16}
- Heterophilic antibodies in a sample have the potential to cause interference in immunoassay systems.^{17, 18} Infrequently, PSA levels may appear elevated due to heterophilic antibodies present in the patient’s serum or plasma or to nonspecific protein binding. If the PSA level is inconsistent with clinical evidence, additional PSA testing is suggested to confirm the result.
- Prostatic massage, ultrasonography and needle biopsy may cause clinically significant elevation in PSA levels. Hormonal therapy may affect PSA expression, therefore a low PSA level after treatment that includes hormonal therapy may not adequately reflect the presence of residual or recurrent disease.
- Some cases of early prostate cancer will not be detected by PSA testing; the same is true for DRE. Biopsy of the prostate is the standard method used to confirm the presence or absence of prostate cancer. For diagnostic purposes, the FastPack® IP Total PSA Immunoassay should always be assessed in conjunction with the patient’s medical history, clinical examination and other findings.

EXPECTED RANGE

Normal, Healthy Male Samples			PSA Level Distribution			
Age (years)	n	Median PSA Level, ng/mL PSA	0 - 4 ng/mL	4.1 - 10.0 ng/mL	10.1 - 30.0 ng/mL	>30.0 ng/mL
30 - 40	9	0.63	9	0	0	0
41 - 50	22	1.24	22	0	0	0
51 - 60	37	0.99	37	0	0	0
61 - 70	25	0.79	24	1	0	0
71 - 80	7	3.42	5	2	0	0
Total male samples	100	0.99	97	3	0	0

Normal, Healthy Male Samples			WHO PSA Level Distribution			
Age (years)	n	Median PSA Level, ng/mL PSA	0 - 3.2 ng/mL	3.3 - 8.0 ng/mL	8.1 - 24.0 ng/mL	>24.0 ng/mL
30 - 40	9	0.50	9	0	0	0
41 - 50	22	0.99	22	0	0	0
51 - 60	37	0.79	37	0	0	0
61 - 70	25	0.63	24	1	0	0
71 - 80	7	2.74	5	2	0	0
Total male samples	100	0.79	97	3	0	0

Normal, Healthy Female Samples			PSA Level Distribution			
	n	Median PSA Level, ng/mL	0 - 4 ng/mL	4.1 - 10.0 ng/mL	10.1 - 30.0 ng/mL	>30.0 ng/mL
	50	0.03	50	0	0	0

Normal, Healthy Female Samples			WHO PSA Level Distribution			
	n	Median PSA Level, ng/mL	0 - 3.2 ng/mL	3.3 - 8.0 ng/mL	8.1 - 24.0 ng/mL	>24.0 ng/mL
	50	0.02	50	0	0	0

Malignant Diseases			PSA Level Distribution			
	n	Median PSA Level, ng/mL	0 - 4 ng/mL	4.1 - 10.0 ng/mL	10.1 - 30.0 ng/mL	>30.0 ng/mL
Kidney	5	0.00	5	0	0	0
Bladder	5	0.44	5	0	0	0
Pancreas	5	0.18	5	0	0	0
Liver	5	0.12	5	0	0	0
Breast	5	0.08	5	0	0	0
Testicles	5	0.61	5	0	0	0
Total Malignant Diseases	30	0.15	30	0	0	0

Malignant Diseases			WHO PSA Level Distribution			
	n	Median PSA Level, ng/mL	0 - 3.2 ng/mL	3.3 - 8.0 ng/mL	8.1 - 24.0 ng/mL	>24.0 ng/mL
Kidney	5	0.00	5	0	0	0
Bladder	5	0.35	5	0	0	0
Pancreas	5	0.14	5	0	0	0
Liver	5	0.12	5	0	0	0
Breast	5	0.10	5	0	0	0
Testicles	5	0.49	5	0	0	0
Total Malignant Diseases	30	0.12	30	0	0	0

Non Malignant Diseases	n	Median PSA Level, ng/mL	PSA Level Distribution			
			0 - 4 ng/mL	4.1 - 10.0 ng/mL	10.1 - 30.0 ng/mL	>30.0 ng/mL
Benign Prostate Hyperplasia	40	1.41	38	2	0	0

Non Malignant Diseases	n	Median PSA Level, ng/mL	WHO PSA Level Distribution			
			0 – 3.2 ng/mL	3.3 - 8.0 ng/mL	8.1 - 24.0 ng/mL	>24.0 ng/mL
Benign Prostate Hyperplasia	40	1.13	38	2	0	0

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The reproducibility of the Total PSA assay was measured by assaying five different levels (n=80, for each level) twenty non-consecutive days using two analyzers and two lots of reagents. The coefficient of variation (%CV) between analyzers, between runs was calculated using analysis of variance.

Lot 1

Sample	Average, ng/mL		Within Run	Between Run	Between Day	Total Precision	%TE	%T'Ea
	Hybritech	WHO	%CV	%CV	%CV	%CV		
1	0.59	0.48	6.739	3.437	0.000	7.565	13.995	33.6
2	2.98	2.38	4.488	1.720	0.000	4.807	18.464	33.6
3	10.91	8.73	7.479	1.015	2.276	7.883	28.001	33.6
4	22.12	17.69	6.888	0.000	4.976	8.498	31.480	33.6
5	31.06	24.85	6.561	4.417	1.664	8.082	19.158	33.6

Lot 2

Sample	Average, ng/mL		Within Run	Between Run	Between Day	Total Precision	%TE	%T'Ea
	Hybritech	WHO	%CV	%CV	%CV	%CV		
1	0.56	0.45	10.123	3.649	1.157	10.823	29.573	33.6
2	2.70	2.16	9.274	2.895	2.952	10.154	22.425	33.6
3	10.36	8.29	9.451	0.000	1.710	9.605	21.788	33.6
4	20.81	16.65	9.597	4.642	7.917	13.278	28.592	33.6
5	31.30	25.04	7.329	3.742	4.296	9.283	22.439	33.6

Spike Recovery:

Serum - Known samples of PSA were added to female serum specimens. The concentration of PSA was determined before and after the addition of the exogenous PSA and the percent recovery was calculated.

Concentration Added (ng/mL)	Observed Concentration (ng/mL)	Recovery (%)
5	4.80	95.9
25	24.64	98.6
40	42.79	107.0

Plasma - Known samples of PSA were added to a normal Heparin pool and a normal EDTA pool. The concentration of PSA was determined and the percentage recoveries were calculated.

Expected ng/mL PSA (Heparin)	Observed ng/mL PSA (Heparin)	Recovery (%)
5.3	5.3	100.0
9.8	9.35	95.4
23.3	29.55	126.8
	Mean	107.4

Expected ng/mL PSA (EDTA)	Observed ng/mL PSA (EDTA)	Recovery (%)
5.0	4.4	88.0
9.5	8.75	92.1
23.0	26.95	117.2
—	Mean	99.1

Dilution Recovery:

Serum - 5 patient pools were diluted 50%, 25% and 12.5% using the zero calibrator and the percentage recoveries were calculated.

Pool #1	Expected Concentration (ng/mL)	Observed Concentration (ng/mL)	Recovery (%)
Neat	0.94	0.94	—
50%	0.47	0.44	94.6
25%	0.23	0.28	120.4
12.5%	0.12	0.13	111.8

Pool #2	Expected Concentration (ng/mL)	Observed Concentration (ng/mL)	Recovery (%)
Neat	57.07	>50	—
50%	28.53	27.91	97.8
25%	14.27	14.6	102.4
12.5%	7.13	7.13	100.0

Pool #3	Expected Concentration (ng/mL)	Observed Concentration (ng/mL)	Recovery (%)
Neat	16.17	16.17	—
50%	8.08	7.62	94.3
25%	4.04	4.43	109.5
12.5%	2.02	2.18	107.9

Pool #4	Expected Concentration (ng/mL)	Observed Concentration (ng/mL)	Recovery (%)
Neat	72.5	>50	—
50%	36.25	35.32	97.4
25%	18.12	17.74	97.9
12.5%	9.06	9.06	100.0

Pool #5	Expected Concentration (ng/mL)	Observed Concentration (ng/mL)	Recovery (%)
Neat	34.79	34.79	—
50%	17.39	18.09	104.0
25%	8.70	8.45	94.21
12.5%	4.35	4.56	104.8

Plasma - Three EDTA and three Heparin patient pools were diluted 50%, 25%, and 12.5% using the appropriate normal pool (i.e., EDTA patient pool diluted with EDTA normal pool) and the percentage recoveries were calculated.

% of Neat EDTA Pool 1	Expected, ng/mL PSA	Observed, ng/mL PSA	Recovery (%)
100	45.5	45.5	—
50	22.9	28.1	122.7
25	11.1	12.2	109.9
12.5	5.7	6.2	108.8
Endogenous Value (ng/mL)	0.31		

% of Neat EDTA Pool 2	Expected, ng/mL PSA	Observed, ng/mL PSA	Recovery (%)
100	13.9	13.9	—
50	7.1	6.9	97.2
25	3.7	3.5	94.6
12.5	2.0	2.0	100.0
Endogenous Value (ng/mL)	0.35		

% of Neat EDTA Pool 3	Expected, ng/mL PSA	Observed, ng/mL PSA	Recovery (%)
100	24.0	24.0	—
50	12.3	11.5	93.5
25	6.4	6.0	93.8
12.5	3.5	3.4	97.1
Endogenous Value (ng/mL)	0.59		

% of Neat Heparin Pool 1	Expected, ng/mL PSA	Observed, ng/mL PSA	Recovery (%)
100	5.3	5.3	—
50	3.1	3.3	106.5
25	2.0	2.1	105.0
12.5	1.4	1.3	92.9
Endogenous Value (ng/mL)	0.85		

% of Neat Heparin Pool 2	Expected, ng/mL PSA	Observed, ng/mL PSA	Recovery (%)
100	6.7	6.7	—
50	3.8	4.5	118.4
25	2.3	2.4	104.3
12.5	1.6	1.5	93.8
Endogenous Value (ng/mL)	0.88		

% of Neat Heparin Pool 3	Expected, ng/mL PSA	Observed, ng/mL PSA	Recovery (%)
100	2.2	2.2	—
50	1.5	1.6	106.7
25	1.2	1.2	100.0
12.5	1.0	1.0	100.0
Endogenous Value (ng/mL)	0.88		

Plasma/Serum Comparison

Clinical samples were used to compare the values obtained from plasma and serum samples from the same patient using the FastPack® IP Total PSA method. The values were evaluated for agreement using Deming regression analysis.

n	Range of Observation (ng/mL)	Intercept (ng/mL)	Slope	r
130	0.0 - 49.1	-0.0597	0.9855	0.97

Method Comparison

189 serum samples were used to compare FastPack® tPSA, 25µL against predicate Beckman Access Hybritech PSA. Evaluation was performed using Passing-Bablok regression analysis.

N	Range of Values (ng/mL)	Intercept (ng/mL)	Slope	R
189	0.06 to 46.8	-0.230	1.082	0.985

N	Range of Values, WHO (ng/mL)	Intercept (ng/mL)	Slope	R
189	0.05 to 37.4	-0.184	1.082	0.985

Equimolarity Studies

Three samples containing approximately 2.5, 5 and 10 ng/mL of PSA and varying proportions of free PSA and PSA ACT complex (0 to 100%) were assayed. Statistical analysis at the 95% confidence interval demonstrated that the proportions of free PSA and PSA-ACT complex did not affect the PSA values (P-value 0.324).

INTERFERING SUBSTANCES

Interfering substances were added to serum pools containing known amounts of PSA. The value obtained for the serum pool with each interfering substance was compared to the value obtained for the serum pool without the interfering substance. These compounds did not show interference at the levels indicated.

Test Compound	Test Concentration	Chemotherapeutic Agents	Concentration
d-Biotin	100 ng/mL	Cyclophosphamide	700 µg/mL
Bilirubin	49 mg/dL	Diethylstilbestrol	1 µg/mL
Hemoglobin	600 mg/dL	Doxorubicin HCl	16 µg/mL
Human IgG	1900 mg/dL	Methotrexate	8 µg/mL
Prostatic Acid Phosphatase (PAP)	1000 ng/dL	Megestrol Acetate	90 µg/mL
Triglycerides	3000 mg/dL	Flutamide	10 µg/mL
		Lupron	100 µg/mL

*Interference was found with Human Serum Albumin at 3 g/dL above the endogenous levels. Interference was found with d-Biotin above 100 ng/mL.

Analytical Sensitivity:

Limit of Blank (LOB), Limit of Detection (LOD), and Limit of Quantitation (LOQ)

The limit of blank (LOB, the highest measurement likely to be observed for a blank sample), limit of detection (LOD, the lowest amount of analyte in a sample that can be detected with type I and II error rates set to 5%), and limit of quantitation (LOQ, the lowest amount of analyte in a sample that can be reliably detected) were determined according to CLSI EP17-A2. In this study the limit of blank was determined from 240 replicate determinations of a blank sample tested on six different FastPack® analyzers using three reagent lots. Raw RLUs from the assays were converted to apparent ng/mL based on the calibration curve for each assay. The LOB was determined as the 228th rank of the sorted distribution of values. This value was 0.01 ng/mL TPSA.

The LOD was estimated from 180 replicate determinations of four low level samples. Per the CLSI EP17-A2 guideline, the parametric LOD calculation was utilized and yielded 0.01 ng/mL TPSA.

The LOQ was determined as the lowest sample which provided <20% CV. The LOQ was set to 0.04 ng/mL.

Analytical Specificity

For the monoclonal antibodies used, the following cross-reactivities were determined: PAP: none; PSA and PSA-ACT are recognized on an equimolar basis.

REFERENCES

- ¹ Lilja H, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines (LMPG): Practice Guidelines and Recommendations for Use of Tumor Markers in the Clinic, Prostate Cancer (Section B), Draft 2006. National Academy of Clinical Biochemistry.
- ² Henttu P, Vihko P. Prostate-specific antigen and human glandular kallikrein: Two kallikreins of the human prostate. *Ann Med* 1994;26:157-164.
- ³ Tewari PC, Bluestein BI. Multiple forms of prostate specific antigen and the influences of immunoassay design on their measurement in patient serum. *J Clin Ligand Assay*, 18 1995;3:186-196.
- ⁴ Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and characterization of different molecular forms of prostatespecific antigen in human seminal fluid. *Clin Chem*. 1995;41(11):1567-1573.
- ⁵ Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. *Clin Lab Invest* 1995;55 Suppl221:32-34.
- ⁶ Scher HI, Kelly WK. Flutamide withdrawal syndrome: Its impact on clinical trials in hormone-refractory prostate cancer. *J Clin Oncol* 1993;11(8):1566-1572.
- ⁷ Partin AW, et al. Prostate specific antigen in the staging of localized prostate cancer: Influence of tumor differentiation, tumor volume and benign hyperplasia. *J Urol* 1990; 143:747-752.
- ⁸ Semjonow A, Brandt B, Oberpenning F, Hertle L. Discrepancies in assays impair the interpretation of prostate-specific antigen. *Urology* 1995;34:303-315.
- ⁹ Armbruster DA. Prostate-Specific Antigen: Biochemistry, Analytical, and Clinical Application [Review]. *Clin Chem* 1993; 39:181.
- ¹⁰ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition; H3-A5; 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS).
- ¹¹ Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- ¹² Primus FJ, Kelly EA, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem*. 1988;34(2):261-264.
- ¹³ Schroff RJ, Foon KA, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-885.
- ¹⁴ Van Duijnhoven HLP, Pérquériaux NCV, van Zon JPHM, Blankenstein MA. Large discrepancy between prostate-specific antigen results from different assays during longitudinal follow-up of a prostate cancer patient. *Clin Chem* 1996;42(4):637-641.
- ¹⁵ Wians FH. The "correct" PSA concentration. *Clin Chem*. 1996;42(11):1882-1885.
- ¹⁶ Cohen RJ, Haffejee Z, Steele GS, Naylar SJ. Advanced prostate cancer with normal serum prostate-specific antigen levels. *Arch Pathol Lab Med* 1994;118:1123-1126.
- ¹⁷ Nahm MH, Hoffman JW. Heteroantibody: Phantom of the immunoassay. *Clin Chem* 1990;36(6):829.
- ¹⁸ Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- ¹⁹ Chybowski FM, Berstralh EJ, Oesterling JE. The effect of digital rectal examination on the serum prostate specific antigen concentration: Results of a randomized study. *J Urol* 1992;148:83-86.
- ²⁰ US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. *Federal Register* 1991; 56(235): 64175-82.

Qualigen, Inc.
Carlsbad, CA 92011 USA
Technical Support
(760) 918-9165
(877) 709-2169

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© 2000 Qualigen, Inc. All rights reserved. Qualigen and FastPack are trademarks or registered trademarks of Qualigen, Inc.
All other trademarks are the property of their respective owners.

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REFERENZEN

- ¹ Lilja H, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines (LMPG): Practice Guidelines and Recommendations for Use of Tumor Markers in the Clinic, Prostate Cancer (Section B), Draft 2006. National Academy of Clinical Biochemistry.
- ² Henttu P, Vihko P. Prostate-specific antigen and human glandular kallikrein: Two kallikreins of the human prostate. *Ann Med* 1994;26:157-164.
- ³ Tewari PC, Bluestein BI. Multiple forms of prostate specific antigen and the influences of immunoassay design on their measurement in patient serum. *J Clin Ligand Assay*, 18 1995;3:186-196.
- ⁴ Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and characterization of different molecular forms of prostatespecific antigen in human seminal fluid. *Clin Chem*. 1995;41(11):1567-1573.
- ⁵ Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. *Clin Lab Invest* 1995;55 Suppl221:32-34.
- ⁶ Scher HI, Kelly WK. Flutamide withdrawal syndrome: Its impact on clinical trials in hormone-refractory prostate cancer. *J Clin Oncol* 1993;11(8):1566-1572.
- ⁷ Partin AW, et al. Prostate specific antigen in the staging of localized prostate cancer: Influence of tumor differentiation, tumor volume and benign hyperplasia. *J Urol* 1990; 143:747-752.
- ⁸ Semjonow A, Brandt B, Oberpenning F, Hertle L. Discrepancies in assays impair the interpretation of prostate-specific antigen. *Urology* 1995;34:303-315.
- ⁹ Armbruster DA. Prostate-Specific Antigen: Biochemistry, Analytical, and Clinical Application [Review]. *Clin Chem* 1993; 39:181.
- ¹⁰ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition; H3-A5; 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS).
- ¹¹ Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- ¹² Primus FJ, Kelly EA, et al. „Sandwich“-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem*. 1988;34(2):261-264.
- ¹³ Schroff RJ, Foon KA, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-885.
- ¹⁴ Van Duijnhoven HLP, Pérquériaux NCV, van Zon JPHM, Blankenstein MA. Large discrepancy between prostate-specific antigen results from different assays during longitudinal follow-up of a prostate cancer patient. *Clin Chem* 1996;42(4):637-641.
- ¹⁵ Wians FH. The „correct“ PSA concentration. *Clin Chem*. 1996;42(11):1882-1885.
- ¹⁶ Cohen RJ, Haffejee Z, Steele GS, Naylar SJ. Advanced prostate cancer with normal serum prostate-specific antigen levels. *Arch Pathol Lab Med* 1994;118:1123-1126.
- ¹⁷ Nahm MH, Hoffman JW. Heteroantibody: Phantom of the immunoassay. *Clin Chem* 1990;36(6):829.
- ¹⁸ Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- ¹⁹ Chybowski FM, Berstralh EJ, Oesterling JE. The effect of digital rectal examination on the serum prostate specific antigen concentration: Results of a randomized study. *J Urol* 1992;148:83-86.
- ²⁰ US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. *Federal Register* 1991; 56(235): 64175-82.

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RÉFÉRENCES

- ¹ Lilja H, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines (LMPG): Practice Guidelines and Recommendations for Use of Tumor Markers in the Clinic, Prostate Cancer (Section B), Draft 2006. National Academy of Clinical Biochemistry.
- ² Henttu P, Vihko P. Prostate-specific antigen and human glandular kallikrein: Two kallikreins of the human prostate. *Ann Med* 1994;26:157-164.
- ³ Tewari PC, Bluestein BI. Multiple forms of prostate specific antigen and the influences of immunoassay design on their measurement in patient serum. *J Clin Ligand Assay*, 18 1995;3:186-196.
- ⁴ Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and characterization of different molecular forms of prostatespecific antigen in human seminal fluid. *Clin Chem*. 1995;41(11):1567-1573.
- ⁵ Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. *Clin Lab Invest* 1995;55 Suppl221:32-34.
- ⁶ Scher HI, Kelly WK. Flutamide withdrawal syndrome: Its impact on clinical trials in hormone-refractory prostate cancer. *J Clin Oncol* 1993;11(8):1566-1572.
- ⁷ Partin AW, et al. Prostate-specific antigen in the staging of localized prostate cancer: Influence of tumor differentiation, tumor volume and benign hyperplasia. *J Urol* 1990; 143:747-752.
- ⁸ Semjonow A, Brandt B, Oberpenning F, Hertle L. Discrepancies in assays impair the interpretation of prostate-specific antigen. *Urology* 1995;34:303-315.
- ⁹ Armbruster DA. Prostate-Specific Antigen: Biochemistry, Analytical, and Clinical Application [Review]. *Clin Chem* 1993; 39:181.
- ¹⁰ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition; H3-A5; 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS).
- ¹¹ Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- ¹² Primus FJ, Kelly EA, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem*. 1988;34(2):261-264.
- ¹³ Schroff RJ, Foon KA, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-885.
- ¹⁴ Van Duijnhoven HLP, Pérquériaux NCV, van Zon JPHM, Blankenstein MA. Large discrepancy between prostate-specific antigen results from different assays during longitudinal follow-up of a prostate cancer patient. *Clin Chem* 1996;42(4):637-641.
- ¹⁵ Wians FH. The "correct" PSA concentration. *Clin Chem*. 1996;42(11):1882-1885.
- ¹⁶ Cohen RJ, Haffejee Z, Steele GS, Naylar SJ. Advanced prostate cancer with normal serum prostate-specific antigen levels. *Arch Pathol Lab Med* 1994;118:1123-1126.
- ¹⁷ Nahm MH, Hoffman JW. Heteroantibody: Phantom of the immunoassay. *Clin Chem* 1990;36(6):829.
- ¹⁸ Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- ¹⁹ Chybowski FM, Berstralh EJ, Oesterling JE. The effect of digital rectal examination on the serum prostate specific antigen concentration: Results of a randomized study. *J Urol* 1992;148:83-86.
- ²⁰ US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. *Federal Register* 1991; 56(235): 64175-82.

Qualigen, Inc.
Carlsbad, CA 92011, États-Unis
Assistance technique:
+1-760-918-9165
+1-877-709-2169

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