

Daily Environmental & QC Log

Take refrigerator and room temperature readings along with relative humidity readings daily. Acceptable ranges for FastPack IP reagent packs and the FastPack IP System analyzer are posted at the beginning of each column. If any of the readings are outside the range, record the corrective action taken (i.e. adjusted temperature control knob and checked testing supplies for deterioration) in the Corrective Action Log.

Month: _____

Day	Refrigerator		Room Temperature		Relative Humidity		Run QC	Clean Equip
	2 – 8° C (35-46° F)		15 – 32° C (59–90° F)		10 – 80%			Procedure Manual
	Temp	Initials	Temp	Initials	%	Initials	Initials	Initials
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								
21								
22								
23								
24								
25								
26								
27								
28								
29								
30								
31								