




Better Patient Experience. Better Patient Retention.

In-office testing in less than 12 minutes allows you to **diagnose and treat in one visit.**

TESTOSTERONE | SHBG | PSA | Vitamin D | TSH | Free T4 | hCG



Why Test In The Office?

THE BENEFITS OF IN-OFFICE TESTING



Patient convenience—no separate lab visit.



Reduces cost.



Results are in hand in less than 12 minutes allowing treatment discussions to take place in person.



Patient visit is more complete and less stressful.



Patient satisfaction equals improved patient retention and long-term success of your practice.

FastPack[®] IP System



Dispense sample



**Place FastPack[®] into
FastPack[®] IP Instrument**




Press START

Qualigen Diagnostics supports you along the way:

- CLIA requirements
- Personnel training
- Individualized financial analysis
- Flexible acquisition options

| ASSAY | SAMPLE TYPE | LOWER LIMIT | UPPER LIMIT | CPT CODE | REIMBURSEMENT |
|--------------|---------------|-------------|-------------|----------|---------------|
| Testosterone | Plasma/Serum | 23 ng/dL | 1600 ng/dL | 84403 | \$25.81 |
| SHBG | Plasma/Serum | 0.80 nmol/L | 174 nmol/L | 84270 | \$21.73 |
| PSA | Plasma/Serum | 0.04 ng/mL | 50 ng/mL | 84153 | \$18.39 |
| Vitamin D | Plasma/Serum | 12.9 ng/mL | 150 ng/mL | 82306 | \$29.60 |
| TSH | Plasma /Serum | 0.13 µIU/mL | 100 µIU/mL | 84443 | \$16.80 |
| free T4 | Serum | 0.4 ng/dL | 6 ng/dL | 84439 | \$9.02 |
| hCG | Plasma/Serum | 1.8 mIU/mL | 1000 mIU/mL | 84702 | \$15.05 |

| SPECIFICATION | DESCRIPTION |
|-------------------------|--|
| Sample Processing | Fully Automatic |
| Sample Material | Plasma/Serum |
| Sample Volume | 100µL per analysis (25µL hCG) |
| Measuring Principle | Chemiluminescent Immunoassay |
| Analysis Throughput | Up to 4 tests/hour |
| Interface | RS-232C |
| Dimensions/Weight | 13 x 9 x 12 inches/29 lbs. |
| Electrical Requirements | 100-240 ± 10% VAC, 1.9 A |
| Operating Conditions | Temperature: 59°F to 90°F Humidity: 10% to 80% Altitude: ≤ 6,562 ft. |



The **FastPack[®] IP System** — Diagnostic technology that gives you the power to provide reliable results in real-time.

QUALIGEN DIAGNOSTICS
2042 Corte Del Nogal, Suite B
Carlsbad, California 92011
877.709.2169 option 3
FastPack@qualigeninc.com

